Hemospray® - Domande frequenti

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- Hookey L, Barkun A, Sultanian R, et al. Successful hemostasis of active lower GI bleeding using a hemostatic powder as monotherapy, combination therapy or rescue therapy. *Gastrointest Endosc.* 2019;89(4):865-871.
- Hussein M, Alzoubaidi D, Weaver M, et al. Use of Hemospray in the treatment of lower gastrointestinal bleeds: outcomes from the international multicenter Hemospray registry. *DDW.* 2020:SU1039
- Hussein M, Alzoubaidi D, O'Donnell M, et al. Hemostatic powder TC-325 treatment of malignancy related upper gastrointestinal bleeds: international registry outcomes. *J Gastroenterol Hepatol.* 2021;36(11):3027-3032.
- Mourad FH, Leong RW. Role of hemostatic powders in the management of lower gastrointestinal bleeding: a review. *J Gastroenterol Hepatol.* 2018;33(8):1445-1453.
- Barkun AN, Adam V, Lu Y, et al. Using Hemospray improves the cost-effectiveness ratio in management of upper gastrointestinal nonvariceal bleeding. *J Clin Gastroenterol.* 2018;52(1):36-44.
- Chen Y, Barkun A. Hemostatic powders in gastrointestinal bleeding: a systematic review. *Gastrointest Endosc Clin N Am.* 2015;25(3):535-552.
- Changela K, Papafragkakis H, Ofori E, et al. Hemostatic powder spray: a new method for managing gastrointestinal bleeding. *Therap Adv Gastroenterol.* 2015;8(3):125-135.
- Sulz MC, Frei R, Meyenberger C, et al. Routine use of Hemospray for gastrointestinal bleeding: prospective two-center experience in Switzerland. *Endoscopy.* 2014;46(7):619-624.
- Yau A, Ou G, Galoport C, et al. Safety and efficacy of Hemospray in upper gastrointestinal bleeding. *Can J Gastroenterol Hepatol.* 2014;28(2):72-76.
- Morris AJ, Smith LA, Stanley A, et al. Hemospray for non-variceal upper gastrointestinal bleeding: results of the Seal Dataset (survey to evaluate the application of Hemospray in the luminal tract). *J Clin Gastroenterol.* E-published December 2013.

Per ulteriori informazioni su Hemospray visitare: [hemospray.cookmedical.com](https://www.hemospray.cookmedical.com)

¹Bustamante-Balén M, Plumé G. Role of hemostatic powders in the endoscopic management of gastrointestinal bleeding. *World J Gastrointest Pathophysiol.* 2014;5(3):284-292. doi:10.4291/wjgp.v5.i3.284.

²Giday SA, Kim Y, Krishnamurthy DM, et al. Long-term randomized controlled trial of a novel nanopowder hemostatic agent (TC-325) for control of severe arterial upper gastrointestinal bleeding in a porcine model. *Endoscopy.* 2011;43(4):296-299. doi:10.1055/s-0030-1256125.

³Alzoubaidi D, Hussein M, Rusu R, et al. Outcomes from an international multicenter registry of patients with acute gastrointestinal bleeding undergoing endoscopic treatment with Hemospray. *Dig Endosc.* 2020;32(1):96-105.

⁴Lau JYW, Pittayanon R, Kwek A, et al. Comparison of a hemostatic powder and standard treatment in the control of active bleeding from upper nonvariceal lesions: a multicenter, noninferiority, randomized trial. *Ann Intern Med.* 2022;175(2):171-178.

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