

It's time to upgrade.

MINC+™

BENCHTOP INCUBATOR



The incubator that IVF clinics
have relied on for two decades
just got better.

We listened to users'
feedback to make
the improvements
that matter to them.



+ Trusted performance

MINC Benchtop Incubators have been an integral part of IVF practice worldwide for 20 years.

+ Humidification

Each chamber is independently humidified.

Gas-flow bubbles and the water level can be observed at the front of the device.



+ DishTrace

The DishTrace platform combines functionality through the incubator touchscreen and PC software to provide a wide range of dish-data management tools.

DishTrace MINC+

DishTrace MINC+ allows you to access and log information via the touchscreen on the center console.

- Add and edit dish records, assign dish locations, and view environmental settings.
- Check in and check out dishes.



DishTrace PC

DishTrace PC software allows you to monitor and manage up to 50 MINC+ incubators.

- View the current status and log the information for up to 50 MINC+ incubators connected to the local PC network.
- Add, edit, search, and view data records for registered dishes in the connected MINC+ incubators.
- Generate and view graphical histories of the dishes, dish locations, and incubator alerts and conditions.
- Generate quality-control reports for all MINC+ activity as PDF or CSV files.
- View dish data that is entered through DishTrace MINC+.



*DishTrace is an aid to logging and viewing the location of dishes within MINC+ incubators. It is not intended to replace laboratory protocols for identifying dishes and their contents.

+ Capacity

Chambers are 50% larger than the previous generation of MINC, with room for 12 Nunc® 4 well dishes.

Media warming chamber holds 6 BD® Falcon® test tubes.

+ Monitoring features

MINC+ continuously monitors critical functions.

Alerts and status icons are visible on the front screens.

+ Peace of mind

Service life is 12 years.

A dedicated service and maintenance team is available.

BD is a registered trademark of Becton, Dickinson and Company. Falcon is a registered trademark of Corning Incorporated. Nunc is a registered trademark of Nunc A/S Corporation.

MINC+™

BENCHTOP INCUBATOR

The MINC+ Benchtop Incubator is intended to be used in the assisted reproduction technology as part of the process to treat infertile patients.

The MINC+ components include:

- MINC+ Benchtop Incubator
- MINC+ user manual
- Two disposable humidification flasks
- Braided Connecting Tube (3 m length)
- Medical grade mains power cord
- USB flash drive with DishTrace PC software
- Ethernet cable
- External alarm connector
- 12 magnetic patient information plates
- Gas inlet T-piece with protective caps
- Test tube cradle
- Front panel USB port cover



Ask your Cook Medical representative about the availability of service and maintenance plans.

Order Number	Reference Part Number	Dimensions mm	Weight kg	Length m
MINC+ Benchtop Incubator				
G44429	K-MINC-2000	556 W x 197 H x 449 D	23.4	–
MINC+ Humidification Flask				
G47317	K-MINC-2000-HF	–	–	–
Braided Connecting Tube				
G53848	K-MINC-BCT-10-300	–	–	3

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Gas supply type: Blend of 6% CO₂, 5% O₂, 89% N₂ (Cook culture system at sea level) or high-purity 6% CO₂ in air (recommended tolerances ± 0.2%)

Gas supply pressure: 150 kPa ± 15 kPa (21.8 psi ± 2.2 psi) (1,500 mbar ± 150 mbar)

Power supply: Universal input 100-240 VAC, 50-60 Hz

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4