Take control.

LithAssist™
SUCTION CONTROL FOR LASER LITHOTRIPSY

www.cookmedical.com

LithAssist™

SUCTION CONTROL FOR LASER LITHOTRIPSY

The LithAssist Suction Control for Laser Lithotripsy is a lightweight, handheld device. It gives you one-handed control of suction during lithotripsy of kidney or bladder stones. Add LithAssist to your procedures and eliminate the need for an assistant who is dedicated to controlling suction.

Use the LithAssist in conjunction with a laser such as the H-30™ Holmium Laser System in order to achieve increased control over all aspects of the procedure. When using the laser to break up the stone, simultaneously suction the fragments and fluids using the LithAssist. This process may ultimately reduce your time in the operating theatre.

The LithAssist has a rigid 11.6 Fr dual-lumen cannula that is attached to an ergonomic handle for control. The smaller 5 Fr lumen provides access for the introduction of a laser fibre for use during lithotripsy of kidney or bladder stones. The larger lumen facilitates the suction of fluid, soft tissue, and stone fragments.

The LithAssist features a trigger mechanism that enables you to manually control the suction strength.

Attach the tubing from an independent vacuum source to the universal adapter at the bottom of the LithAssist.

Benefits

- The LithAssist gives you single-handed control of the suction. This allows you to suction fragments and fluids simultaneously during laser lithotripsy without the help of an assistant.
- It works with your existing holmium laser system. This eliminates the need to make an additional capital investment for percutaneous nephrolithotomy (PCNL).
- It locks the laser fibre in place. This allows you to maneuver the tip of the fibre with precision and control.
- It avoids the need to make direct contact with the stone during the procedure. This may reduce trauma to the urothelial tissue.

1. The ergonomic handle with a locking trigger controls suction.
2. The stainless steel cannula delivers the laser fibre to the stone and facilitates the suction of fluids and fragments.
3. The Tuohy-Borst adapter secures the laser fibre in place.
4. The Luer lock is for syringe introduction to flush out clots and clogs.

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