



# Spectrum<sup>®</sup>

MINOCYCLINE/RIFAMPIN CATHETER



Proven lowest catheter-related  
bloodstream infection rates.



# Contents

---

## CATHETER SETS

7	4.0 Fr Double Lumen
8	5.0 Fr Double Lumen/5.0 Fr Triple Lumen
9	7.0 Fr Triple Lumen
10	8.0 Fr Double Lumen
11	10.0 Fr Five Lumen

---

## FULL SPECTRUM TRAYS

15	4.0 Fr Double Lumen
16	5.0 Fr Double Lumen/5.0 Fr Triple Lumen
17	7.0 Fr Triple Lumen
18	8.0 Fr Double Lumen
19	9.0 Fr Triple Lumen
20	10.0 Fr Five Lumen

---

## IHI TRAYS

23	4.0 Fr Double Lumen
24	5.0 Fr Double Lumen/5.0 Fr Triple Lumen
25	7.0 Fr Triple Lumen
26	8.0 Fr Double Lumen
27	9.0 Fr Triple Lumen
28	10.0 Fr Five Lumen

---

## 29-30 CUSTOM TRAY PROGRAM

---

# PROCESS + SPECTRUM CATHETERS

To combat CRBSIs, we impregnate our polyurethane catheters with the clinically proven combination of minocycline and rifampin. More than 21 peer-reviewed articles and meta-analyses recognize Spectrum as superior—it's evidence of protection no other process or technology can match. Furthermore, over 10 years of clinical use and a study of more than 500,000 catheter days show no evidence of bacterial resistance.<sup>1</sup> Cook Spectrum® catheters give your patients the best protection possible against CRBSI and may also provide significant economic benefits for your hospital.



### THE RIGHT COMBINATION

Minocycline+rifampin is the most rigorously studied and synergistic combination of antibiotics available on a catheter for reducing CRBSIs through two distinct mechanisms of action. Unlike most antibiotics, this combination has the ability to penetrate the biofilm that forms on all indwelling catheters.<sup>2</sup>



### UNMATCHED EVIDENCE

Two decades of evidence, including more than 21 peer-reviewed studies and four meta-analyses, confirm minocycline+rifampin catheters are the most effective tool available to provide broad-spectrum protection against gram-positive, gram-negative and fungal infections<sup>3,4</sup> in both short- and long-term use. It's evidence of protection no other process or technology can match.



### NO EVIDENCE OF RESISTANCE

More than 10 years of clinical use has shown no evidence that M+R catheters lead to bacterial resistance, and a seven-year study of over 500,000 catheter days confirms these results.<sup>1</sup> A separate study indicates that a facility using M+R catheters may have shown a decrease in vancomycin usage.<sup>6</sup>



### WORTH SWITCHING

In a challenging clinical environment, a hospital that switches to M+R catheters may expect to see a decrease in CRBSI rates, attributable mortality and CRBSI-related costs. Even high-performing hospitals can switch to Spectrum to drive incremental improvement in CRBSI rates and may still achieve substantial reductions in mortality and expenses.







Set includes:

- Power-injectable\* Cook Spectrum central venous catheter with hydrophilic coating and Luer lock end caps\*
- Safe-T-J® double flexible-tipped wire guide with cm markings\*
- Dilator
- Introducer needle†
- FEP catheter introducer needle
- Movable suture wing‡
- Disposable 6 mL syringe
- Needleless injection caps
- Red sharps needle holder

\* Available only on 7, 8, 9 and 10 Fr catheters.

† EchoTip® echogenic needle is standard with 4 Fr double-lumen (standard material), 5 Fr double-, 5 Fr triple- and 8 Fr double-lumen CVCs.

‡ Selected catheter lengths only.

PROCESS + SPECTRUM

Proven Lowest  
Infection Rates



## 4.0 Fr Double Lumen




Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Soft Material</b>						
G10041	C-SDLM-401J-PED-ABRM	4.0	18	5	.018	21 TW
G10039	C-SDLM-401J-ABRM	4.0	18	8	.018	21 TW
G10438	C-SDLM-401J-RSC-ABRM	4.0	18	12	.018	21 TW
G10452	C-SDLM-401J-LSC-ABRM	4.0	18	15	.018	21 TW
<b>Standard Material</b>						
G43984	C-UDLM-401J-PED-ABRM-HC	4.0	18	5	.018	21 TW
G43983	C-UDLM-401J-ABRM-HC	4.0	18	8	.018	21 TW
G43985	C-UDLM-401J-RSC-ABRM-HC	4.0	18	12	.018	21 TW
G43986	C-UDLM-401J-LSC-ABRM-HC	4.0	18	15	.018	21 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



## 5.0 Fr Double Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
G43996	C-UDLM-501J-PED-ABRM-HC	5.0	16	5	.018	21 TW
G43995	C-UDLM-501J-ABRM-HC	5.0	16	8	.018	21 TW
G43997	C-UDLM-501J-RSC-ABRM-HC	5.0	16	12	.018	21 TW
G43998	C-UDLM-501J-LSC-ABRM-HC	5.0	16	15	.018	21 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## 5.0 Fr Triple Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
G44006	C-UTLM-501J-PED-ABRM-HC	5.0	16	5	.025	19 TW
G44005	C-UTLM-501J-ABRM-HC	5.0	16	8	.025	19 TW
G44007	C-UTLM-501J-RSC-ABRM-HC	5.0	16	12	.025	19 TW
G44008	C-UTLM-501J-LSC-ABRM-HC	5.0	16	15	.025	19 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## 7.0 Fr Triple Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Regular</b>						
G47825	C-UTLM-701J-ABRM-HC-RD	7.0	13	15	.032	18 TW
G47829	C-UTLM-701J-RSC-ABRM-HC-RD	7.0	13	20	.032	18 TW
G47827	C-UTLM-701J-LSC-ABRM-HC-RD	7.0	13	25	.032	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## 8.0 Fr Double Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Regular</b>						
G50817	C-UDLM-801J-ABRM-HC	8.0	12	15	.032	18 TW
G50818	C-UDLM-801J-RSC-ABRM-HC	8.0	12	20	.032	18 TW
<b>Minipuncture™</b>						
G44063	C-UDLM-801J-MNP-ABRM-HC	8.0	12	15	.025	20 XTW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.




## 10.0 Fr Five Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
G51537	C-UQLM-1001J-ABRM-HC-RD	10.0	10	15	.032	18 TW
G51538	C-UQLM-1001J-LSC-ABRM-HC-RD	10.0	10	25	.032	18 TW
G51539	C-UQLM-1001J-RSC-ABRM-HC-RD	10.0	10	20	.032	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



Proven Lowest  
Infection Rates



# Full Spectrum Trays



ChloroPrep is a registered trademark of CareFusion 2200, Inc.  
BakSnap is a registered trademark of DuoProSS Meditech Corporation.



### Full Spectrum Tray contains:

- Power-injectable\* Cook Spectrum® central venous catheter with hydrophilic coating and Luer lock end caps\*
- Safe-T-J® double flexible-tipped wire guide with cm markings\*
- EchoTip® echogenic introducer needles, 5 cm and 7 cm
- Catheter insertion checklist
- BakSnap® safety syringes
- 10 mL prefilled sodium chloride syringes
- 25 gage needle
- 22 gage needle
- FEP catheter introducer needle
- Transducer tubing
- Dilator
- Lidocaine and lidocaine label
- Filter straw
- Tinted ChloraPrep® One-Step
- Full-body fenestrated drape with clear window
- Gauze sponges
- Disposable safety scalpel
- Straight needle with suture or curved needle with suture and needle holder
- Disposable syringes
- Needleless injection caps
- Locking sharps needle holder
- Movable suture wing\*\*

PROCESS + SPECTRUM

Proven Lowest  
Infection Rates



\*Available only on 7, 8, 9 and 10 Fr catheters.

\*\*Selected catheter lengths only.

## 4.0 Fr Double Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Soft Material – Curved Needle with Suture and Holder</b>						
G29418	C-SDLMY-401J-PED-ABRM-FST	4.0	18	5	.018	21 TW
G29413	C-SDLMY-401J-ABRM-FST	4.0	18	8	.018	21 TW
G29420	C-SDLMY-401J-RSC-ABRM-FST	4.0	18	12	.018	21 TW
G29416	C-SDLMY-401J-LSC-ABRM-FST	4.0	18	15	.018	21 TW
<b>Standard Material – Curved Needle with Suture and Holder</b>						
G43992	C-UDLMY-401J-PED-ABRM-HC-FST	4.0	18	5	.018	21 TW
G43991	C-UDLMY-401J-ABRM-HC-FST	4.0	18	8	.018	21 TW
G43993	C-UDLMY-401J-RSC-ABRM-HC-FST	4.0	18	12	.018	21 TW
G43994	C-UDLMY-401J-LSC-ABRM-HC-FST	4.0	18	15	.018	21 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



## 5.0 Fr Double Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Curved Needle with Suture and Holder</b>						
G43982	C-UDLMY-501J-PED-ABRM-HC-FST	5.0	16	5	.018	21 TW
G43981	C-UDLMY-501J-ABRM-HC-FST	5.0	16	8	.018	21 TW
G44003	C-UDLMY-501J-RSC-ABRM-HC-FST	5.0	16	12	.018	21 TW
G44004	C-UDLMY-501J-LSC-ABRM-HC-FST	5.0	16	15	.018	21 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## 5.0 Fr Triple Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Curved Needle with Suture and Holder</b>						
G44014	C-UTLMY-501J-PED-ABRM-HC-FST	5.0	16	5	.025	19 TW
G44013	C-UTLMY-501J-ABRM-HC-FST	5.0	16	8	.025	19 TW
G44015	C-UTLMY-501J-RSC-ABRM-HC-FST	5.0	16	12	.025	19 TW
G44016	C-UTLMY-501J-LSC-ABRM-HC-FST	5.0	16	15	.025	19 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



## 7.0 Fr Triple Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Curved Needle with Suture and Holder</b>						
G44126	C-UTLMY-701J-ABRM-HC-FST-A-RD	7.0	13	15	.032	18 TW
G44128	C-UTLMY-701J-RSC-ABRM-HC-FST-A-RD	7.0	13	20	.032	18 TW
G44127	C-UTLMY-701J-LSC-ABRM-HC-FST-A-RD	7.0	13	25	.032	18 TW
<b>Straight Needle with Suture</b>						
G47831	C-UTLMY-701J-ABRM-HC-FST-RD	7.0	13	15	.032	18 TW
G47833	C-UTLMY-701J-RSC-ABRM-HC-FST-RD	7.0	13	20	.032	18 TW
G47832	C-UTLMY-701J-LSC-ABRM-HC-FST-RD	7.0	13	25	.032	18 TW
<b>Curved Needle with Suture, Holder and Minipuncture™ Access</b>						
G47835	C-UTLMY-701J-MNP-ABRM-HC-FST-A-RD	7.0	13	15	.025	20 TW
G47839	C-UTLMY-701J-RSC-MNP-ABRM-HC-FST-A-RD	7.0	13	20	.025	20 TW
G47837	C-UTLMY-701J-LSC-MNP-ABRM-HC-FST-A-RD	7.0	13	25	.025	20 TW
<b>Pigtail Catheter – Curved Needle with Suture and Holder</b>						
G47845	C-UTLMY-701J-190-PIG-ABRM-FST-A-RD	7.0	13	15	.035	18 TW
G47846	C-UTLMY-701J-190-PIG-RSC-ABRM-FST-A-RD	7.0	13	20	.035	18 TW
<b>Pigtail Catheter – Straight Needle with Suture</b>						
G47844	C-UTLMY-701J-190-PIG-ABRM-FST-RD	7.0	13	15	.035	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## 8.0 Fr Double Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Curved Needle with Suture, Holder and Minipuncture Access</b>						
G44071	C-UDLMY-801J-LSC-MNP-ABRM-HC-FST-A	8.0	12	25	.025	20 XTW
<b>Straight Needle with Suture and Minipuncture Access</b>						
G44067	C-UDLMY-801J-RSC-MNP-ABRM-HC-FST	8.0	12	20	.025	20 XTW
G44068	C-UDLMY-801J-LSC-MNP-ABRM-HC-FST	8.0	12	25	.025	20 XTW
<b>Curved Needle with Suture and Holder</b>						
G50819	C-UDLMY-801J-ABRM-HC-FST-A	8.0	12	15	.032	18 TW
G50820	C-UDLMY-801J-RSC-ABRM-HC-FST-A	8.0	12	20	.032	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



## 9.0 Fr Triple Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Curved Needle with Suture and Holder</b>						
G51586	C-UTLMY-901J-ABRM-HC-FST-A-RD	9.0	11	15	.032	18 TW
G51587	C-UTLMY-901J-RSC-ABRM-HC-FST-A-RD	9.0	11	20	.032	18 TW
G51588	C-UTLMY-901J-LSC-ABRM-HC-FST-A-RD	9.0	11	25	.032	18 TW
<b>Straight Needle with Suture</b>						
G51585	C-UTLMY-901J-ABRM-HC-FST-RD	9.0	11	15	.032	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## 10.0 Fr Five Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Curved Needle with Suture and Holder</b>						
G51543	C-UQLMY-1001J-ABRM-HC-FST-A-RD	10.0	10	15	.032	18 TW
G51545	C-UQLMY-1001J-RSC-ABRM-HC-FST-A-RD	10.0	10	20	.032	18 TW
G51544	C-UQLMY-1001J-LSC-ABRM-HC-FST-A-RD	10.0	10	25	.032	18 TW
<b>Straight Needle with Suture</b>						
G51540	C-UQLMY-1001J-ABRM-HC-FST-RD	10.0	10	15	.032	18 TW
G51542	C-UQLMY-1001J-RSC-ABRM-HC-FST-RD	10.0	10	20	.032	18 TW
G51541	C-UQLMY-1001J-LSC-ABRM-HC-FST-RD	10.0	10	25	.032	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



# Maximal Sterile Barrier Precautions Trays

The Spectrum procedural tray has been uniquely designed to support your hospital's process bundle. It contains several components that comply with national recommendations for patient safety from AHRQ, CDC, IHI, Joint Commission, OSHA and SHEA Compendium.





#### Tray contains:

- Power-injectable\* Cook Spectrum® central venous catheter with hydrophilic coating and Luer lock end caps\*
- Safe-T-J® double flexible-tipped wire guide with cm markings\*
- EchoTip® echogenic introducer needles, 5 cm and 7 cm
- Catheter insertion checklist
- BakSnap® safety syringes
- 25 gage needle
- 22 gage needle
- FEP catheter introducer needle
- Transducer tubing
- Dilator
- Lidocaine and lidocaine label
- Filter straw
- Tinted Chloraprep® One-Step
- Full-body fenestrated drape with clear window
- Gauze sponges
- Disposable safety scalpel
- Straight needle with suture or curved needle with suture and needle holder
- Disposable syringes
- Needleless injection caps
- Prefilled sodium chloride syringes, 10 mL
- Locking sharps needle holder
- Movable suture wing\*\*
- Cap, mask and gown

PROCESS + SPECTRUM

Proven Lowest  
Infection Rates



\*Available only on 7, 8, 9 and 10 Fr catheters.

\*\*Selected catheter lengths only.

## 4.0 Fr Double Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Soft Material – Straight Needle with Suture and Holder</b>						
G48049	C-SDLMY-401J-PED-ABRM-IHI-FST	4.0	18	5	.018	21 TW
G48047	C-SDLMY-401J-ABRM-IHI-FST	4.0	18	8	.018	21 TW
G48050	C-SDLMY-401J-RSC-ABRM-IHI-FST	4.0	18	12	.018	21 TW
G48048	C-SDLMY-401J-LSC-ABRM-IHI-FST	4.0	18	15	.018	21 TW
<b>Standard Material – Straight Needle with Suture</b>						
G48059	C-UDLMY-401J-PED-ABRM-HC-IHI-FST	4.0	18	5	.018	21 TW
G48053	C-UDLMY-401J-ABRM-HC-IHI-FST	4.0	18	8	.018	21 TW
G48061	C-UDLMY-401J-RSC-ABRM-HC-IHI-FST	4.0	18	12	.018	21 TW
G48056	C-UDLMY-401J-LSC-ABRM-HC-IHI-FST	4.0	18	15	.018	21 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



## 5.0 Fr Double Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Straight Needle with Suture and Holder</b>						
G48069	C-UDLMY-501J-PED-ABRM-HC-IHI-FST	5.0	16	5	.018	21 TW
G48064	C-UDLMY-501J-ABRM-HC-IHI-FST	5.0	16	8	.018	21 TW
G48076	C-UDLMY-501J-RSC-ABRM-HC-IHI-FST	5.0	16	12	.018	21 TW
G48067	C-UDLMY-501J-LSC-ABRM-HC-IHI-FST	5.0	16	15	.018	21 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## 5.0 Fr Triple Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Straight Needle with Suture</b>						
G48122	C-UTLMY-501J-PED-ABRM-HC-IHI-FST	5.0	16	5	.025	19 TW
G48113	C-UTLMY-501J-ABRM-HC-IHI-FST	5.0	16	8	.025	19 TW
G48124	C-UTLMY-501J-RSC-ABRM-HC-IHI-FST	5.0	16	12	.025	19 TW
G48117	C-UTLMY-501J-LSC-ABRM-HC-IHI-FST	5.0	16	15	.025	19 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



## 7.0 Fr Triple Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Curved Needle with Suture and Holder</b>						
G49803	C-UTLMY-701J-ABRM-HC-IHI-FST-A-RD	7.0	13	15	.032	18 TW
G49804	C-UTLMY-701J-RSC-ABRM-HC-IHI-FST-A-RD	7.0	13	20	.032	18 TW
G49807	C-UTLMY-701J-LSC-ABRM-HC-IHI-FST-A-RD	7.0	13	25	.032	18 TW
<b>Straight Needle with Suture</b>						
G49802	C-UTLMY-701J-ABRM-HC-IHI-FST-RD	7.0	13	15	.032	18 TW
G49805	C-UTLMY-701J-RSC-ABRM-HC-IHI-FST-RD	7.0	13	20	.032	18 TW
G49806	C-UTLMY-701J-LSC-ABRM-HC-IHI-FST-RD	7.0	13	25	.032	18 TW
<b>Curved Needle with Suture, Holder and Minipuncture™ Access</b>						
G49809	C-UTLMY-701J-MNP-ABRM-HC-IHI-FST-A-RD	7.0	13	15	.025	20 XTW
G49811	C-UTLMY-701J-RSC-MNP-ABRM-HC-IHI-FST-A-RD	7.0	13	20	.025	20 XTW
G49813	C-UTLMY-701J-LSC-MNP-ABRM-HC-IHI-FST-A-RD	7.0	13	25	.025	20 XTW
<b>Straight Needle with Suture and Minipuncture Access</b>						
G49808	C-UTLMY-701J-MNP-ABRM-HC-IHI-FST-RD	7.0	13	15	.025	20 XTW
G49827	C-UTLMY-701J-RSC-MNP-ABRM-HC-IHI-FST-RD	7.0	13	20	.025	20 XTW
G49812	C-UTLMY-701J-LSC-MNP-ABRM-HC-IHI-FST-RD	7.0	13	25	.025	20 XTW
<b>Pigtail Catheter – Curved Needle with Suture and Holder</b>						
G49815	C-UTLMY-701J-190-PIG-ABRM-IHI-FST-A-RD	7.0	13	15	.035	18 TW
G49817	C-UTLMY-701J-190-PIG-RSC-ABRM-IHI-FST-A-RD	7.0	13	20	.035	18 TW
G49819	C-UTLMY-701J-190-PIG-LSC-ABRM-IHI-FST-A-RD	7.0	13	25	.035	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## 8.0 Fr Double Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Curved Needle with Suture and Holder</b>						
G50821	C-UDLMY-801J-ABRM-HC-IHI-FST-A	8.0	12	15	.032	18 TW
G50822	C-UDLMY-801J-RSC-ABRM-HC-IHI-FST-A	8.0	12	20	.032	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



## 9.0 Fr Triple Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Curved Needle with Suture and Holder</b>						
G51589	C-UTLMY-901J-ABRM-HC-IHI-FST-A-RD	9.0	11	15	.032	18 TW
G51590	C-UTLMY-901J-RSC-ABRM-HC-IHI-FST-A-RD	9.0	11	20	.032	18 TW
G51591	C-UTLMY-901J-LSC-ABRM-HC-IHI-FST-A-RD	9.0	11	25	.032	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## 10.0 Fr Five Lumen



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>Curved Needle with Suture and Holder</b>						
G51546	C-UQLMY-1001J-ABRM-HC-IHI-FST-A-RD	10.0	10	15	.032	18 TW
G51548	C-UQLMY-1001J-RSC-ABRM-HC-IHI-FST-A-RD	10.0	10	20	.032	18 TW
G51547	C-UQLMY-1001J-LSC-ABRM-HC-IHI-FST-A-RD	10.0	10	25	.032	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



## Enhance your process bundle with customized protection.

### Volume Commitment

If your hospital orders at least 75 Spectrum trays per month from Cook Medical, it qualifies for the program.

### How to Order

Orders will be placed via an electronic checklist. Please contact your local Cook representative to complete the checklist.

### Timeline

You will receive your first order approximately 10 weeks after the electronic checklist is completed.

### Cardinal Presource<sup>®</sup> Partnership

If the custom tray program doesn't meet your needs, Cook Medical now partners with Cardinal Health's . . . Presource<sup>®</sup> to provide Central Line Insertion Packs.



## Spectrum Custom Tray Program

### Standard Components

- One multi-lumen central venous catheter\*
- 18 gage EchoTip® introducer needles, 5 and 7 cm lengths
- 18 gage catheter access percutaneous entry needle, 7 cm length
- For Minipuncture™: 20 gage EchoTip and 20 gage catheter access percutaneous entry needle, 5 cm length
- Transducer tube
- 6 mL introducer syringe
- Ten 4 x 4 inch gauze pads
- Two 2 x 2 inch gauze pads
- Locking sharps container
- #11 blade safety scalpel
- Patient sticker, labels, IFU



### Custom Tray Options

#### Wire Guide (Choose 1)

- .032 inch marked wire guide (not available with pigtail catheters)
- .035 inch PTFE coated – all sizes
- .025 inch Minipuncture wire guide, 7 Fr only

#### Maximal Sterile Barrier (Choose Quantity of 0, 1 or 2)

- Cap, gown and mask with eyeshield (all or none)
- Drape (choose quantity of 0 or 1)
  - 68 x 72 inch drape
  - 72 x 122 inch drape

#### Skin Prep (Choose Quantity of 0, 1 or 2)

- 3 mL ChloraPrep® One-Step, orange tint

#### Flush Syringes (Choose Yes or No)

- 10 mL prefilled sodium chloride syringes (quantity varies by number of lumens)

#### Local Anesthetic (Choose Quantity of 0 or 1 for Each)

- 5 mL 1% lidocaine with label
- 22 gage needle
- 25 gage needle
- Filter straw

#### Syringes (Choose Quantity of 0, 1 or 2)

- 3 mL Luer lock syringe
- 3 mL BakSnap® safety syringe

#### Suture and Needle (Choose 1)

- 3-0 braided silk with straight needle
- 3-0 braided silk with curved needle and driver
- 3-0 monofilament nylon with straight needle
- 3-0 monofilament nylon with curved needle and driver

#### Dressings (Choose Quantity of 0 or 1 for Each)

- 1 inch diameter BioPatch®
- 3.5 x 4.5 inch Tegaderm® IV

#### Injection Caps (Choose 0 or Number of Lumens for Each)

- CLC2000® positive pressure
- MicroCLAVE® neutral displacement

#### Insertion Checklist (Choose Yes or No)

- Standard Cook Medical-generated checklist

\*Spectrum catheters available in 4, 5, 7, 8, 9 and 10 Fr sizes only.  
Uncoated catheters available in 4, 5, 7, 9 and 10 Fr sizes only.

## References

1. Ramos ER, Reitzel R, Jiang Y, et al. Clinical effectiveness and risk of emerging resistance associated with prolonged use of antibiotic-impregnated catheters: more than 0.5 million catheter days and 7 years of clinical experience. *Crit Care Med.* 2011;39(1):1-7.
2. Raad I, Hanna H, Jiang Y, et al. Comparative activities of daptomycin, linezolid, and tigecycline against catheter-related methicillin-resistant *Staphylococcus bacteremic* isolates embedded in biofilm. *Antimicrob Agents Chemother.* 2007;51(5):1656-1660.
3. Schuerer DJ, Mazuski JE, Buchman TG, et al. Catheter-related bloodstream infection rates in minocycline/rifampin vs. chlorhexidine/silver sulfadiazine-impregnated central venous catheters - results of a 46 month study. *Crit Care Med.* 2008;36(12)(suppl):A199-A208. Abstract 454.
4. Raad I, Darouiche R, Hachem R, et al. The broad-spectrum activity and efficacy of catheters coated with minocycline and rifampin. *J Infect Dis.* 1996;173(2):418-424.
5. Brooks K, Dauenhauer S, Nelson M. Comparison of an untreated vs. silver/chlorhexidine vs. rifampin/minocycline central venous catheter in reducing catheter-related bloodstream infections. Abstract presented at: APIC 28th Annual Educational Conference and International Meeting; June 10-14, 2001; Seattle, WA.

## Customer Service Centers

AORTIC INTERVENTION  
CRITICAL CARE  
ENDOSCOPY  
INTERVENTIONAL RADIOLOGY  
LEAD MANAGEMENT  
PERIPHERAL INTERVENTION  
SURGERY  
UROLOGY  
WOMEN'S HEALTH

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

