



## COMPONENT ORDERING INFORMATION



Order Number	Reference Part Number	Distal Body Diameter mm	Common Iliac Segment Length mm	External Iliac Segment Length mm	Introducer Sheath ID Fr
<b>Iliac Bifurcation Grafts</b>					
G38612	ZBIS-10-45-41	10	45	41	20
G38614	ZBIS-10-61-41	10	61	41	20
G38613	ZBIS-10-45-58	10	45	58	20
G38615	ZBIS-10-61-58	10	61	58	20
G38616	ZBIS-12-45-41	12	45	41	20
G38618	ZBIS-12-61-41	12	61	41	20
G38617	ZBIS-12-45-58	12	45	58	20
G38344	ZBIS-12-61-58	12	61	58	20



Order Number	Reference Part Number	Distal Iliac Leg Diameter mm	Iliac Leg Working Length mm	Introducer Sheath ID Fr
<b>Zenith® Spiral-Z™ Iliac Legs</b>				
G55236	ZSLE-16-39-ZT	16	39	14
G55237	ZSLE-16-56-ZT	16	56	14
G55238	ZSLE-16-74-ZT	16	74	14
G55239	ZSLE-16-90-ZT	16	90	14

AORTIC INTERVENTION  
CRITICAL CARE  
ENDOSCOPY  
INTERVENTIONAL RADIOLOGY  
LEAD MANAGEMENT  
PERIPHERAL INTERVENTION  
SURGERY  
UROLOGY  
WOMEN'S HEALTH



### Customer Service Centers

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
Poland: +48 223060159, [pl.orders@cookmedical.com](mailto:pl.orders@cookmedical.com)  
Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

