



Zenith Dissection

Graft Specifications and Components

Zenith[®] Dissection


ENDOVASCULAR SYSTEM

Zenith® Dissection


ENDOVASCULAR SYSTEM

GRAFT SPECIFICATIONS


Proximal component




Proximal tapered component



Proximal double tapered component








Z-Trak Plus® Introduction System

22-34 mm grafts: 20 Fr (6.7 mm) ID/7.7 mm OD

36-42 mm grafts: 22 Fr (7.3 mm) ID/8.5 mm OD



Note: All measurements are in millimeters. The black numbers on the illustrations indicate stent lengths. The red numbers indicate gap lengths. Gap lengths are approximate.



Zenith® Dissection

ENDOVASCULAR SYSTEM

PROXIMAL COMPONENTS

Proximal component

Order Number	Reference Part Number	Component Type	Aortic Vessel OD mm	Diameter mm	Length mm	Introducer Sheath ID/OD Fr (mm)/mm
22 mm Diameter						
G47437	ZDEG-P-22-79-PF-US	P	20	22	79	20 (6.7)/7.7
G47436	ZDEG-P-22-117-PF-US	P	20	22	117	20 (6.7)/7.7
24 mm Diameter						
G47439	ZDEG-P-24-79-PF-US	P	21	24	79	20 (6.7)/7.7
G47438	ZDEG-P-24-117-PF-US	P	21	24	117	20 (6.7)/7.7
26 mm Diameter						
G47441	ZDEG-P-26-79-PF-US	P	22-23	26	79	20 (6.7)/7.7
G47440	ZDEG-P-26-136-PF-US	P	22-23	26	136	20 (6.7)/7.7
28 mm Diameter						
G47444	ZDEG-P-28-82-PF-US	P	24	28	82	20 (6.7)/7.7
G47442	ZDEG-P-28-142-PF-US	P	24	28	142	20 (6.7)/7.7
G47443	ZDEG-P-28-202-PF-US	P	24	28	202	20 (6.7)/7.7
30 mm Diameter						
G47447	ZDEG-P-30-82-PF-US	P	25-27	30	82	20 (6.7)/7.7
G47445	ZDEG-P-30-142-PF-US	P	25-27	30	142	20 (6.7)/7.7
G47446	ZDEG-P-30-202-PF-US	P	25-27	30	202	20 (6.7)/7.7
32 mm Diameter						
G47450	ZDEG-P-32-82-PF-US	P	28-29	32	82	20 (6.7)/7.7
G47448	ZDEG-P-32-142-PF-US	P	28-29	32	142	20 (6.7)/7.7
G47449	ZDEG-P-32-202-PF-US	P	28-29	32	202	20 (6.7)/7.7
34 mm Diameter						
G47453	ZDEG-P-34-79-PF-US	P	30	34	79	20 (6.7)/7.7
G47451	ZDEG-P-34-154-PF-US	P	30	34	154	20 (6.7)/7.7
G47452	ZDEG-P-34-204-PF-US	P	30	34	204	20 (6.7)/7.7
36 mm Diameter						
G47456	ZDEG-P-36-79-PF-US	P	31-32	36	79	22 (7.3)/8.5
G47454	ZDEG-P-36-154-PF-US	P	31-32	36	154	22 (7.3)/8.5
G47455	ZDEG-P-36-204-PF-US	P	31-32	36	204	22 (7.3)/8.5
38 mm Diameter						
G47459	ZDEG-P-38-79-PF-US	P	33-34	38	79	22 (7.3)/8.5
G47457	ZDEG-P-38-154-PF-US	P	33-34	38	154	22 (7.3)/8.5
G47458	ZDEG-P-38-204-PF-US	P	33-34	38	204	22 (7.3)/8.5
40 mm Diameter						
G47462	ZDEG-P-40-83-PF-US	P	35-36	40	83	22 (7.3)/8.5
G47460	ZDEG-P-40-164-PF-US	P	35-36	40	164	22 (7.3)/8.5
G47461	ZDEG-P-40-218-PF-US	P	35-36	40	218	22 (7.3)/8.5
42 mm Diameter						
G47465	ZDEG-P-42-83-PF-US	P	37-38	42	83	22 (7.3)/8.5
G47463	ZDEG-P-42-164-PF-US	P	37-38	42	164	22 (7.3)/8.5
G47464	ZDEG-P-42-218-PF-US	P	37-38	42	218	22 (7.3)/8.5

PROXIMAL SINGLE / DOUBLE TAPERED COMPONENTS

Proximal tapered component

Order Number	Reference Part Number	Component Type	Aortic Vessel OD mm	Proximal Diameter mm	Distal Diameter mm	Length mm	Introducer Sheath ID/OD Fr (mm)/mm
32/28 mm Diameter							
G47466	ZDEG-PT-32-162-PF-US	PT	28-29	32	28	162	20(6.7)/7.7
G47467	ZDEG-PT-32-202-PF-US	PT	28-29	32	28	202	20(6.7)/7.7
34/30 mm Diameter							
G47470	ZDEG-PT-34-159-PF-US	PT	30	34	30	159	20(6.7)/7.7
G47471	ZDEG-PT-34-199-PF-US	PT	30	34	30	199	20(6.7)/7.7
36/32 mm Diameter							
G47474	ZDEG-PT-36-159-PF-US	PT	31-32	36	32	159	22(7.3)/8.5
G47475	ZDEG-PT-36-199-PF-US	PT	31-32	36	32	199	22(7.3)/8.5
38/34 mm Diameter							
G47478	ZDEG-PT-38-154-PF-US	PT	33-34	38	34	154	22(7.3)/8.5
G47479	ZDEG-PT-38-204-PF-US	PT	33-34	38	34	204	22(7.3)/8.5
40/36 mm Diameter							
G47482	ZDEG-PT-40-160-PF-US	PT	35-36	40	36	160	22(7.3)/8.5
G47483	ZDEG-PT-40-210-PF-US	PT	35-36	40	36	210	22(7.3)/8.5
42/38 mm Diameter							
G47486	ZDEG-PT-42-160-PF-US	PT	37-38	42	38	160	22(7.3)/8.5
G47487	ZDEG-PT-42-210-PF-US	PT	37-38	42	38	210	22(7.3)/8.5

Proximal double tapered component

Order Number	Reference Part Number	Component Type	Aortic Vessel OD mm	Proximal Diameter mm	Distal Diameter mm	Length mm	Introducer Sheath ID/OD Fr (mm)/mm
32/24 mm Diameter							
G47468	ZDEG-PT-32-24-158-PF-US	PT	28-29	32	24	158	20 (6.7)/7.7
G47469	ZDEG-PT-32-24-196-PF-US	PT	28-29	32	24	196	20 (6.7)/7.7
34/26 mm Diameter							
G47472	ZDEG-PT-34-26-156-PF-US	PT	30	34	26	156	20 (6.7)/7.7
G47473	ZDEG-PT-34-26-194-PF-US	PT	30	34	26	194	20 (6.7)/7.7
36/28 mm Diameter							
G47476	ZDEG-PT-36-28-159-PF-US	PT	31-32	36	28	159	22 (7.3)/8.5
G47477	ZDEG-PT-36-28-199-PF-US	PT	31-32	36	28	199	22 (7.3)/8.5
38/30 mm Diameter							
G47480	ZDEG-PT-38-30-159-PF-US	PT	33-34	38	30	159	22 (7.3)/8.5
G47481	ZDEG-PT-38-30-199-PF-US	PT	33-34	38	30	199	22 (7.3)/8.5
40/32 mm Diameter							
G47485	ZDEG-PT-40-32-165-PF-US	PT	35-36	40	32	165	22 (7.3)/8.5
G47484	ZDEG-PT-40-32-205-PF-US	PT	35-36	40	32	205	22 (7.3)/8.5
42/34 mm Diameter							
G47488	ZDEG-PT-42-34-160-PF-US	PT	37-38	42	34	160	22 (7.3)/8.5
G47489	ZDEG-PT-42-34-210-PF-US	PT	37-38	42	34	210	22 (7.3)/8.5

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.


Zenith® Dissection

ENDOVASCULAR SYSTEM

ENDOVASCULAR STENTS


Order Number	Reference Part Number	Aortic Vessel OD mm	Diameter mm	Length mm	Stents	Introducer Sheath ID/OD Fr (mm)/mm
36 mm Diameter						
G47491	ZDES-36-80-US	20-30	36	80	4	16 (5.3)/6.0
G47492	ZDES-36-120-US	20-30	36	120	6	16 (5.3)/6.0
G47493	ZDES-36-180-US	20-30	36	180	9	16 (5.3)/6.0
46 mm Diameter						
G47494	ZDES-46-80-US	31-38	46	80	4	16 (5.3)/6.0
G47495	ZDES-46-120-US	31-38	46	120	6	16 (5.3)/6.0
G47496	ZDES-46-185-US	31-38	46	185	9	16 (5.3)/6.0

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Z-Trak Plus Introduction System for dissection stents

16 Fr (5.3 mm) ID/6.0 mm OD



INDICATIONS FOR USE

	The Zenith Dissection Endovascular System (Zenith TX2 Dissection Endovascular Graft with Pro-Form and Zenith Dissection Endovascular Stent) is indicated for the endovascular treatment of patients with Type B aortic dissection. The Zenith TX2 Dissection Endovascular Graft with Pro-Form is intended to seal the entry tear(s) and to exclude aneurysms associated with chronic dissections. The Zenith Dissection Endovascular Stent is intended to be used as a distal component to provide support to delaminated segments of non-aneurysmal aorta with dissection distal to a Zenith TX2 Dissection Endovascular Graft with Pro-Form.			
Indications for use				
Length	Dissection stent lengths 80, 120, 180, 185 used in combination with the sealing stent graft TX2 Proximal Component can potentially allow treatment of the entire descending thoracic aorta and the abdominal aorta.			
Aortic fixation site diameter	Dissection grafts: 20–38 mm (measured outer-wall-to-outer-wall) Dissection stents: 20–38 mm (measured outer-wall-to-outer-wall)			
Angulation/curvature	• Radius of curvature > 35 mm along the length of aorta that is intended to be treated (dissection graft) • Localized angulation < 45 degrees (dissection graft). • Radius of curvature > 35 mm along the length of aorta that is intended to be treated (dissection stent) • Localized angulation < 45 degrees (dissection stent)			
Access	Adequate femoral/iliac access compatible with the required introduction systems			
Introducer diameters	22-34 mm grafts: ID 20 Fr (6.7 mm) OD 7.7 mm	36-42 mm grafts: ID 22 Fr (7.3 mm) OD 8.5 mm	36, 46 mm stents: ID 16 Fr (5.3 mm) OD 6.0 mm	
Contraindications	• Patients with known sensitivities or allergies to stainless steel, polyester, nitinol, polypropylene, or gold • Patients with a systemic infection who may be at increased risk of endovascular stent/graft infection			



Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, is.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com

AI-ESC-IR-OHNS-PI-RH-SUR-A4