



Zenith® Iliac Branch

COMPONENT ORDERING INFORMATION



Order Number	Reference Part Number	Distal Iliac Leg Diameter mm	Working Length of Iliac Leg mm	Introducer Sheath ID Fr
Zenith® Spiral-Z® Iliac Leg Grafts				
G55236	ZSLE-16-39-ZT	16	39	14
G55237	ZSLE-16-56-ZT	16	56	14
G55238	ZSLE-16-74-ZT	16	74	14
G55239	ZSLE-16-90-ZT	16	90	14



Order Number	Reference Part Number	Distal Body Diameter mm	Common Iliac Segment Length mm	External Iliac Segment Length mm	Introducer Sheath ID Fr
Zenith® Iliac Branch Grafts					
G60340	ZBIS-10-45-41-US	10	45	41	20
G60342	ZBIS-10-61-41-US	10	61	41	20
G60344	ZBIS-12-45-41-US	12	45	41	20
G60345	ZBIS-12-61-41-US	12	61	41	20
G60341	ZBIS-10-45-58-US*	10	45	58	20
G60343	ZBIS-10-61-58-US*	10	61	58	20
G60346	ZBIS-12-45-58-US*	12	45	58	20
G60347	ZBIS-12-61-58-US*	12	61	58	20

*Nonstock

Getinge iCast[®] covered stent system specifications

Cook Order Number	Cook RPN/ Getinge Iliac Code	Diameter (mm)	Length (mm)	Post Nominal Pressure 8 ATM mm		Post Burst Rate Pressure 12 ATM mm		Introducer Sheath Fr
				OD	Length	OD	Length	
iCast covered stent–80 cm Catheter Length								
G60492	48838	8	38	8.1	34.7	8.5	34.7	7
G60493	48859	8	59	8.0	56.0	8.4	56.5	7
G60494	48938	9	38	8.9	33.7	9.3	32.7	7
G60495	48959	9	59	8.9	54.6	9.3	54.0	7
G60513	48038	10	38	10.0	30.8	10.4	30.9	7*
iCast covered stent–120 cm Catheter Length								
G60509	42838	8	38	8.1	34.7	8.5	34.7	7
G60510	42859	8	59	8.0	56.0	8.4	56.5	7
G60511	42938	9	38	8.9	33.7	9.3	32.7	7
G60512	42959	9	59	8.9	54.6	9.3	54.0	7
G60513	42038	10	38	10.0	30.8	10.4	30.9	7*

*If using a 10 mm device with a Cook introducer sheath, it is recommended to use an 8 Fr sheath to minimize potential for resistance.

Protected by the following international and/or U.S. patent(s): <http://patents.getinge.com>. CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician. Getinge and **GETINGE** * are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates in the United States or other countries. iCast is a registered trademark of Atrium Medical Corporation. Copyright 2024 Atrium Medical Corporation. All rights reserved. Refer to Getinge's Instructions for Use for current indications, warnings, contraindications, and precautions.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4