


# Access and navigate the ureter.




## Amplatz

FIXED CORE WIRE GUIDE

### Stiffness options

The Amplatz wire guides are available in stiff, extra stiff, and ultra stiff options. The stiffness of each type of wire guide depends on variations in its mandril, coil, and inner safety wire.



## Amplatz

FIXED CORE WIRE GUIDE

The Amplatz Fixed Core Wire Guides are used to gain ureteral access and establish a tract. Use an Amplatz wire guide to place, replace, or exchange medical devices during urological procedures. The wire guides are PTFE coated and are available in three stiffness options as well as flexible and double-flexible tip configurations.

These wire guides are not intended for PTCA use.

Order Number	Reference Part Number	Diameter inch	Length cm	Tip Configuration	Coating
<b>Stiff</b>					
G53211	AS-035145	.035	145	5 cm flexible	PTFE
G53556	AS-038145	.038	145	5 cm flexible	PTFE
G53557	AS-038145-DF	.038	145	5 cm double flexible	PTFE
<b>Extra Stiff</b>					
G53558	AES-025145	.025	145	5 cm flexible	PTFE
G53559	AES-035145	.035	145	5 cm flexible	PTFE
G53560	AES-038145	.038	145	5 cm flexible	PTFE
G53563	AES-035145-DF	.035	145	5 cm double flexible	PTFE
G53564	AES-038145-DF	.038	145	5 cm double flexible	PTFE
<b>Ultra Stiff</b>					
G53566	AUS-035145	.035	145	5 cm flexible	PTFE
G53567	AUS-038145	.038	145	5 cm flexible	PTFE
G53568	AUS-035145-DF	.035	145	5 cm double flexible	PTFE
G53569	AUS-038145-DF	.038	145	5 cm double flexible	PTFE

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

### Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Poland: +48 223060159, [pl.orders@cookmedical.com](mailto:pl.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)

EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)

#### Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

#### Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



©CC, URO-A4