
Shiley® Flexible Adult Tracheostomy Tube



Features and benefits

- The beveled tip is designed to facilitate percutaneous insertion.
- The TaperGuard™ cuff reduces fluid leakage by 96% at a lower cuff pressure.*
- The tube is designed to conform to a patient's anatomy with a transparent, soft flange, which can promote airflow around the stoma.*

Customer service

For questions, or to **reorder** the individual tracheostomy tubes, please contact **Medtronic's customer service** at **800.635.5267.**

*Shiley adult flexible tracheostomy tube with TaperGuard cuff, disposable inner cannula claims test report. 10122514, revision C (2014).

The Blue Rhino® G2-Multi Percutaneous Tracheostomy Introducer Sets and Trays have been designed to create an optimal fit with the Shiley Flexible Adult Tracheostomy Tubes.

Recommended loading dilator sizes and associated tracheostomy tubes


Loading dilator	Compliant ISO-standard tracheostomy tube sizes
6.5 mm	6.5 mm
7.0 mm	7.0 mm
7.5 mm	7.5 mm
8.0 mm	8.0 mm
8.5 mm	8.5 mm
9.0 mm	9.0 mm
10 mm	10 mm

Compatibility testing was performed with Shiley Flex and Shiley Evac tracheostomy tubes. When using another appropriately sized tracheostomy tube, ensure that the tracheostomy tube's tip fits snugly on the dilator.

Shiley and TaperGuard are registered trademarks of Covidien LP.



Shiley[®] Flexible Evac Adult Tracheostomy Tube



Features and benefits

- The tube is designed to enable the removal of secretions above the cuff.
- The beveled tip is designed to facilitate percutaneous insertion.
- The TaperGuard™ cuff reduces fluid leakage by 96% at a lower cuff pressure.*
- The tube is designed to conform to a patient's anatomy with a transparent, soft flange, which can promote airflow around the stoma.*

Customer service

For questions, or to **reorder** the individual tracheostomy tubes, please contact **Medtronic's customer service** at **800.635.5267.**

*Shiley adult flexible tracheostomy tube with TaperGuard cuff, disposable inner cannula claims test report. 10122514, revision C (2014).

The Blue Rhino® G2-Multi Percutaneous Tracheostomy Introducer Sets and Trays have been designed to create an optimal fit with the Shiley Flexible Adult Tracheostomy Tubes.

Recommended loading dilator sizes and associated tracheostomy tubes

Loading dilator	Compliant ISO-standard tracheostomy tube sizes
6.5 mm	6.5 mm
7.0 mm	7.0 mm
7.5 mm	7.5 mm
8.0 mm	8.0 mm
8.5 mm	8.5 mm
9.0 mm	9.0 mm
10 mm	10 mm

Compatibility testing was performed with Shiley Flex and Shiley Evac tracheostomy tubes. When using another appropriately sized tracheostomy tube, ensure that the tracheostomy tube's tip fits snugly on the dilator.

