




# Precise placement for one-lung ventilation



ENDOBONCHIAL  
BLOCKER SETS



## Arndt

### ENDOBRONCHIAL BLOCKER


The Arndt Endobronchial Blocker Set is intended to differentially intubate a patient's bronchus in order to isolate the left or right lung for procedures that require one-lung ventilation. The 5 French Arndt Endobronchial Blocker is indicated for paediatric populations in children one year and older. All other sizes are for adult use only.

- The adjustable guide loop is coupled with an appropriately sized bronchoscope to enable precise placement.
- The low-pressure, high-volume balloon creates surface area contact with the inner bronchial wall while minimising potential trauma to the bronchus.
- The distance markings indicate how far the blocker has been advanced.
- The suction adapter included with the set offers an additional method of lung deflation.

| Order Number                            | Reference Part Number | Catheter Fr | Catheter Length cm | Recommended ETT* Diameter mm | Average Balloon Inflation Volume mL | Balloon Configuration |
|---|-----------------------|-------------|--------------------|------------------------------|-------------------------------------|-----------------------|
| <b>Arndt Endobronchial Blocker Sets</b> |                       |             |                    |                              |                                     |                       |
| G44109                                  | C-AEBS-5.0-50-SPH-AS  | 5           | 50                 | ≥ 4.5                        | 0.5 - 2.0                           | Paediatric spherical  |
| G44114                                  | C-AEBS-7.0-65-SPH-AS  | 7           | 65                 | ≥ 6.0                        | 2.0 - 6.0                           | Adult spherical       |
| G44120                                  | C-AEBS-9.0-78-SPH-AS  | 9           | 78                 | ≥ 7.5                        | 4.0 - 8.0                           | Adult spherical       |

\*Endotracheal tube

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details. Please see product risk information in the IFU at [cookmedical.eu](http://cookmedical.eu).



## Cohen

### ENDOBRONCHIAL BLOCKER

The Cohen Endobronchial Blocker Set is intended to differentially intubate a patient's bronchus in order to isolate the left or right lung for procedures that require one-lung ventilation. The product is intended for use by physicians trained and experienced in one-lung ventilation. Standard techniques for placement of endobronchial blockers should be employed.

- The unique deflecting tip provides a wide range of motion to precisely direct the blocker under bronchoscopic guidance.
- The low-pressure, high-volume balloon creates surface area contact with the inner bronchial wall while minimising potential trauma to the bronchus.
- Distance markings provide visual guidance by indicating how far the blocker has been advanced.
- The suction adapter included with the set offers an additional method of lung deflation.

| Order Number                           | Reference Part Number | Catheter Fr | Catheter Length cm | Recommended ETT* Diameter mm | Average Balloon Inflation Volume mL | Balloon Configuration |
|--|-----------------------|-------------|--------------------|------------------------------|-------------------------------------|-----------------------|
| <b>Cohen Endobronchial Blocker Set</b> |                       |             |                    |                              |                                     |                       |
| G44122                                 | C-AEBTS-9.0-65-SPH-AS | 9           | 65                 | ≥ 7.0                        | 6.0 - 9.0                           | Adult spherical       |

\*Endotracheal tube

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details. Please see product risk information in the IFU at [cookmedical.eu](http://cookmedical.eu).

## Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)

Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
Poland: +48 223060159, [pl.orders@cookmedical.com](mailto:pl.orders@cookmedical.com)  
Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)

EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)

### Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: [customersupport@cookmedical.com](mailto:customersupport@cookmedical.com)

### Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



CC, URO-A4