

From cytology to histology: Advance your EBUS diagnostic procedures.



EchoTip ProCore[®]

ENDOBONCHIAL HD ULTRASOUND BIOPSY NEEDLE


At Cook, we were among the first to recognise the clinical impact of Endobronchial ultrasound (EBUS) procedures. That's why we focused our expertise as market leaders in ultrasound needles to create a full range of EchoTip ProCore EBUS needles to meet your needs. Now you can potentially obtain cytology and histology in your EBUS procedures, even in the challenging pulmonary hilar region. Our goal is simple: to give you the histological information to help you deliver the best patient care.




Olympus



Pentax



25 gage



22 gage

EchoTip ProCore®

ENDOBONCHIAL HD ULTRASOUND BIOPSY NEEDLE

For use with Olympus: Used with an ultrasound bronchoscope for fine needle biopsy (FNB) of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree or gastrointestinal tract.

For use with Pentax: Used with an ultrasound bronchoscope for fine needle biopsy (FNB) of submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree.

This device is supplied sterile and is disposable. This single-use device is not designed for re-use.

- The needle features a core trap designed for receiving tissue in the delicate pulmonary area.
- It has high-definition dimpling pattern for needle visibility.
- The naturally contoured handle offers needle control and stability.
- The sheath length is adjustable.
- On scopes with metal Luer hubs, the removable sheath adjuster allows direct connection to the scope.
- On scopes with metal non-Luer hubs, the sturdy scope adapter design allows:
 - Secure connection and procedural stability
 - Suction and injection capability due to the airtight valve

Order Number	Reference Part Number	Needle gauge	Sheath Fr	Adjustable Needle Extension cm	Stylet Tip	Minimum Accessory Channel Diameter mm
G34281	ECHO-HD-22-EBUS-O-C*	22	4.1	0-5	recessed ball tip	2
G34279	ECHO-HD-22-EBUS-P-C**	22	4.1	0-5	recessed ball tip	2
G34282	ECHO-HD-25-EBUS-O-C*	25	4.1	0-5	bevelled	2
G34280	ECHO-HD-25-EBUS-P-C**	25	4.1	0-5	bevelled	2

*Needle attaches to scope adapter on Olympus Endobronchial Ultrasound scope with non-Luer metal hub.

**Needle attaches to Luer metal hub on Pentax Endobronchial Ultrasound scope.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.

Customer Service

EU website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, emea.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 21612723770, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com
 EDI: cookmedical.com/edi.do

Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 Email: customersupport@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 Email: cau.custserv@cookmedical.com



CC, URO-A4