

Enhance your diagnostic capabilities.



EchoTip Ultra[®]

ENDOBONCHIAL HD ULTRASOUND NEEDLE

Cook's EchoTip Ultra Endobronchial HD Ultrasound Needle is designed to enhance your diagnostic and staging capabilities. The high echogenicity of the EchoTip Ultra aids in needle visualisation, ensuring the needle tip is within the target to facilitate sample acquisition. The ergonomic handle and scope adapter can also optimise procedural performance, stability, and security.



Olympus



Pentax

EchoTip Ultra®

ENDOBONCHIAL HD ULTRASOUND NEEDLE

For use with Olympus: This device is used to sample targeted submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree or gastrointestinal tract through the accessory channel of an ultrasound bronchoscope for fine needle aspiration (FNA).

For use with Pentax: This device is used to sample targeted submucosal and extramural lesions or tissue within or adjacent to the tracheobronchial tree through the accessory channel of an ultrasound bronchoscope for fine needle aspiration (FNA).

This device is supplied sterile and is disposable. This single-use device is not designed for re-use.

- The needle has a high-definition dimpling pattern for needle visibility.
- The naturally contoured handle offers needle control and stability.
- On scopes with metal Luer hubs, the removable sheath adjuster allows direct connection to the scope.
- On scopes with metal non-Luer hubs, the sturdy scope adaptor design allows:
 - Secure connection and procedural stability
 - Suction and injection capability due to the airtight valve
- An integrated sheath adjuster is included for compatibility with a full range of ultrasound bronchoscopes.

Order Number	Reference Part Number	Needle gauge	Sheath Fr	Adjustable Needle Extension cm	Stylet Tip	Minimum Accessory Channel Diameter mm
G52011	ECHO-HD-22-EBUS-O*	22	4.1	0-5	bevelled	2
G53408	ECHO-HD-22-EBUS-P**	22	4.1	0-5	bevelled	2
G24062	ECHO-HD-25-EBUS-O*	25	4.1	0-5	bevelled	2
G24061	ECHO-HD-25-EBUS-P**	25	4.1	0-5	bevelled	2

*Needle attaches to scope adapter on Olympus Endobronchial Ultrasound scope with metal non-Luer hub.

**Needle attaches to metal Luer hub on Pentax Endobronchial Ultrasound scope.

Please follow any step-by-step instructions carefully, including the risk, technical information, Device Name and intended use provided in the IFU at cookmedical.eu. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

Customer Service

EU website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, emea.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 21612723770, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com
 EDI: cookmedical.com/edi.do

Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 Email: customersupport@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 Email: cau.custserv@cookmedical.com



CC, URO-A4