

# Take your EBUS staging procedures to the next level.

**EchoTip ProCore®**  
Endobronchial HD  
Ultrasound Biopsy Needle

**EchoTip Ultra®**  
Endobronchial HD  
Ultrasound Needle



Order Number	Reference Part Number	Needle gauge	Needle Extension cm	Sheath Fr	Stylet Tip
G52011	ECHO-HD-22-EBUS-O	22	0-5	4.1	beveled
G53408	ECHO-HD-22-EBUS-P	22	0-5	4.1	beveled
G24062	ECHO-HD-25-EBUS-O	25	0-5	4.1	beveled
G24061	ECHO-HD-25-EBUS-P	25	0-5	4.1	beveled
G34281	ECHO-HD-22-EBUS-O-C	22	0-5	4.1	recessed ball tip
G34279	ECHO-HD-22-EBUS-P-C	22	0-5	4.1	recessed ball tip
G34282	ECHO-HD-25-EBUS-O-C	25	0-5	4.1	beveled
G34280	ECHO-HD-25-EBUS-P-C	25	0-5	4.1	beveled

*Minimum accessory channel diameter 2.0 mm. 5 and 10 mL vacuum locking syringe included.*

*G52011, G24062 G34281, and G34282 attach to scope adapter on Olympus Endobronchial Ultrasound scope with non-Luer metal hub.*

*G53408, G24061, G34279, and G34280 attach to Luer metal hub on Pentax Endobronchial Ultrasound scope.*

Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

WITH Please follow any step-by-step instructions carefully, including the risk and technical information provided in the IFU. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

