

# Band with confidence.



## 4,6,10 Shooter™

SAEED MULTI-BAND LIGATOR

After leading the market for more than a decade, the reliable Saeed Multi-Band Ligator is now available with Universal-fit barrels that streamline your order information and simplify inventory control. And, by using bands that demonstrate lower slippage rates when compared to other bands on the market<sup>1</sup>, the Saeed 4, 6, 10 Shooter can potentially reduce the risk of rebleeding. We've taken reliability and streamlined it.

<sup>1</sup>"In-vitro bench studies have shown differences in the resting inner diameters and radial retraction force of bands on different ligation devices. The larger resting inner diameter and a significantly lower retraction force in neoprene bands, compared to rubber bands predisposed these bands to early slippage." J. M. Polski, E. M. Brunt, Z. A. Saeed (2001). Chronology of Histological Changes after Band Ligation of Esophageal Varices in Humans. *Endoscopy*, 33, 443-447.

<sup>2</sup>Ibid, p.446



## 4,6,10 Shooter™

SAEED MULTI-BAND LIGATOR

Used to endoscopically ligate esophageal varices at or above the gastroesophageal junction or to ligate internal hemorrhoids.

Supplied non-sterile and is disposable - intended for single use only.

- Universal-fit barrels streamline ligation procedures and inventory control
- Opti-Vu® barrel features two-string trigger cord and recessed bands – providing maximum endoscopic view while maintaining suction
- Clear band signals when only one band remains to be deployed
- Wide array of barrel sizes to meet a full range of procedural needs

*"The larger resting inner diameter and a significantly lower retraction force in neoprene bands, compared to rubber bands predisposed these bands to early slippage."*<sup>2</sup>

Order Number	Reference Part Number	Endoscope Outer Diameter mm	Number of Bands	Trigger Cord Length cm
G31915	MBL-U-4*	8.6 - 11.3	4	122
G31917	MBL-U-6*	8.6 - 11.3	6	122
G31919	MBL-U-10*	8.6 - 11.3	10	122
G31918	MBL-U-6-F*	8.6 - 11.3	6	142
G24654	MBL-6-OV*	9.5 - 11.5	6	122
G22613	MBL-4	9.5 - 13	4	122
G22553	MBL-6	9.5 - 13	6	122
G22614	MBL-10	9.5 - 11.5	10	122

Order Number	Reference Part Number	Endoscope Outer Diameter mm	Number of Bands	Trigger Cord Length cm
G25122	MBL-4-XL	11 - 14	4	122
G25119	MBL-6-XL	11 - 14	6	122
G25221	MBL-6-XL-C	11 - 14	6	190
G22692	MBL-6-F	9.5 - 13	6	142

\*Kit includes Opti-Vu barrel w/ ligation bands and trigger cord, loading catheter, precision control handle and irrigation adapter. Minimum accessory channel 2.8 mm.

**Caution:** Product contains natural rubber latex, which may cause an allergic reaction.

**Note:** The universal-fit barrels are not compatible with XL scopes.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

### Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)

EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)

#### Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

#### Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: [cau.custserv@cookmedical.com-A4](mailto:cau.custserv@cookmedical.com-A4)



AI, ESC, IR, OHNS, PI, RH, SUR-A4