


Choose intuitive products
that were designed with
patients in mind.



Entuit[®] Gastrostomy BR


BALLOON RETENTION FEEDING TUBE

Entuit Gastrostomy BR


Place the feeding tube with the slim, smooth, in-line balloon shaft that is comfortable for your patient.

The **Entuit Gastrostomy BR Balloon Retention Feeding Tube** sizes match perfectly with the initial placement gastrostomy sets to streamline placement procedures. (Available both with and without ENFit® Connection*.)

The ENFit Connection* helps enhance patient safety by reducing misconnections.



The inner lumen maintains the same diameter even when the balloon is inflated.




The tubes are clearly marked with fill volume and Fr size.

The snug bolster stays in place.

The balloon material is designed for durability at appropriate fill volumes.

The silicone material is soft for patient comfort.



The unique, integrated design of the in-line balloon shaft on 14 Fr and larger tubes reduces the need for overdistension during placement and replacement procedures.




*Also available without ENFit Connection.

Intended to provide gastric access for enteral feeding, medication administration, and decompression through an established gastrointestinal stoma tract. Indicated for use in percutaneous placement of an enteral feeding tube in adult and paediatric patients that require enteral feeding, medication administration, or decompression through an established gastrointestinal stoma tract.


Entuit® Gastrostomy BR

BALLOON RETENTION FEEDING TUBE
WITH ENFIT® CONNECTION



Entuit® Gastrostomy BR

BALLOON RETENTION FEEDING TUBE



Order Number	Reference Part Number	Fr	Marked Length cm	Balloon Volume mL	Used with Entuit Start Initial Placement Gastrostomy Set* sizes
Entuit Gastrostomy BR Balloon Retention Feeding Tube with ENFit Connection					
G44780	SBRD-12-ENF	12	10	5-7	small
G44781	SBRD-14-ENF	14	10	5-7	small
G44782	SBRD-16-ENF	16	10	5-7	medium
G44783	SBRD-18-ENF	18	10	5-10	medium
G44784	SBRD-20-ENF	20	10	20	large
G44785	SBRD-22-ENF	22	10	20	large
G44786	SBRD-24-ENF	24	10	20	large
Entuit Gastrostomy BR Balloon Retention Feeding Tube					
G35614	SBRD-12	12	10	5-7	small
G35615	SBRD-14	14	10	5-7	small
G35616	SBRD-16	16	10	5-7	medium
G35617	SBRD-18	18	10	5-10	medium
G35618	SBRD-20	20	10	20	large
G35619	SBRD-22	22	10	20	large
G35620	SBRD-24	24	10	20	large

*Must be ordered separately.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

ENFit is a federally registered trademark of GEDSA in multiple jurisdictions throughout the world. Used with permission.



ENFit®

GEDSA (The Global Enteral Device Supplier Association) is a non-profit trade association established in part to help introduce the new ISO standard connector and facilitate adoption of the new ENFit connectors within the healthcare community. Learn more at stayconnected.org.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com-A4



AI, ESC, IR, OHNS, PI, RH, SUR-A4