


Mark the spot
for radiation therapy.



EchoTip[®] Ultra
FIDUCIAL NEEDLE



Image courtesy of Marc Giovannini, Endoscopic Unit,
Paoli-Calmettes Institute, Marseille, France.

Delivering targeted therapy



We realized that both endoscopists and radiation oncologists are always looking for better ways to more accurately and efficiently deliver targeted therapy to their patients. That's why we developed the EchoTip Ultra Fiducial Needle. Now, in one procedure, you can place up to four fiducial markers that can be utilized for future therapeutic procedures.

Fiducial markers potentially:

- Increase the accuracy of radiation treatment
- Allow controlled, targeted radiation dose administration which may reduce radiation exposure for the patient
- Provide for improved patient outcomes because you can target the precise area
- Preserve healthy tissue by only targeting diseased tissue




With the EchoTip Ultra Fiducial Needle:

- Place up to four fiducial markers, potentially saving valuable procedure time when compared to traditional methods of manually loading fiducials
- Controlled, precise placement can be achieved due to the tactile feel when the marker is deployed



Place multiple fiducials in one procedure

The EchoTip Ultra Fiducial Needle is the first preloaded EUS-guided fiducial needle, so you do not have to pause to load fiducials after each placement. Based on benchtop modeling¹ the flexible sheath design of the EchoTip Ultra Fiducial needle allows implicit placement of up to four preloaded fiducials in challenging anatomical locations under endoscopic ultrasound. This can potentially enhance your procedural efficiency, which can benefit you, radiation oncologists and your patients.



¹Data on file at Cook Medical.

Product Specifications

Order Number	Reference Part Number	Needle Gage	Sheath Size Fr	Fiducial Length mm	Fiducial Width mm	Number of Fiducials per Device	Minimum Accessory Channel mm
G03793	ECHO-22-F	22	5.2	5	.43	4	3.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Collaborative education and training

We are always seeking to provide stimulating and productive educational opportunities for clinicians to enhance their knowledge of new technologies and learn about exciting new procedural techniques. Cook Medical's Vista programs encourage physician and industry collaboration to improve patient care by sharing best practices. Vista programs foster collaboration among the best and the brightest gastroenterologists in the world who all share one common goal: Improving patient care.

We also share our experience in EUS to facilitate training and education through online resources in the following locations:

 youtube.com/cookmedicalendoscopy

Endoscopic Ultrasound
Video Library - cookmedical.com

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com-A4



AI, ESC, IR, OHNS, PI, RH, SUR-A4