


Cannulate with confidence.



Tracer Metro[®] Direct[™]

WIRE GUIDE

With a kink-resistant core and highly radiopaque tip, the Tracer Metro Direct wire guide provides the ideal balance of maneuverability and visibility you need to confidently cannulate and navigate even difficult ductal anatomies.



Tracer Metro[®] Direct™

WIRE GUIDE

Used to assist in cannulation and to aid in bridging difficult strictures during ERCP.

Placement of the Tracer Metro Direct allows for continual access to the desired duct and facilitates exchange of compatible wire guided accessories.

- Vividly colored spirals confirm wire guide position and movement
- 5 cm hydrophilic coated radiopaque tip provides smooth navigation with enhanced radiopacity
- Kink-resistant core delivers superior maneuverability
- Color Connection facilitates easy identification of compatible accessories

Order Number	Reference Part Number	Wire Guide Length cm	Wire Guide Diameter inch	Radiopaque Tip	Color Connection
G55705	METII-21-260	260	.021	Straight	Black
G26862	METII-21-480	480	.021	Straight	Black
G31913	METII-25-260	260	.025	Straight	Blue
G25785	METII-25-480	480	.025	Straight	Blue
G31914	METII-35-260	260	.035	Straight	Purple
G25686	METII-35-480	480	.035	Straight	Purple
G25687	METII-35-480-A	480	.035	Angled	Purple

May be left in place during sphincterotomy. **Caution:** Not for use with metal tip devices.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com-A4



AI, ESC, IR, OHNS, PI, RH, SUR-A4