Access, orientation and security.



Conveniently fitted with a DomeTip[®] and dedicated wire guide and injection lumens, the Tri-Tome pc delivers the exceptional access and injection capability you need to allow you to focus on better patient outcomes. And, on our Protector[™] model, we've added insulation on the proximal portion of Tri-Tome's cutting wire to reduce the incidence of contact between the cutting wire and the endoscope.





Tri-Tome pc[®]

Used for endoscopic cannulation of the ductal system and for sphincterotomy.

- DomeTip is designed with anatomy in mind for smooth, easy access
- Custom 3-D forming wire for orientation
- Protector[™] insulation on proximal portion of cutting wire decreases incidence of cutting wire contact with scope, reducing the risk of damage to the cutting wire and scope (available on TRI-25M-P)

Order Number	Reference Part Number	Tip	Cutting Wire Length mm	Cutting Wire Type	Wire Guide* Diameter inch	Color Connection
Tri-Tome pc [®] and Tri-Tome pc [®] Protector™						
G22583	TRI-20	DomeTip	20	Braided	.035	Purple
G22584	TRI-20M	DomeTip	20	Monofilament	.035	Purple
G22554	TRI-25	DomeTip	25	Braided	.035	Purple
G22555	TRI-25M	DomeTip	25	Monofilament	.035	Purple
G26826	TRI-25M-SLT	DomeTip	25	Monofilament	.035	Purple
G22556	TRI-30	DomeTip	30	Braided	.035	Purple
G22557	TRI-30M	DomeTip	30	Monofilament	.035	Purple
G49197	TRI-25M-P	DomeTip	25	Monofilament	.035	Purple

Minimum accessory channel 2.8 mm. 7 FR catheter. *Recommend Acrobat* 2 Calibrated Tip Wire Guide, sold separately. Active cord available separately. Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu EDI: cookmedical.eu/edi Distributors: +353 61239240, ssc.distributors@cookmedical.com Austria: +43 179567121, oe.orders@cookmedical.com Belgium: +32 27001633, be.orders@cookmedical.com Denmark: +45 38487607, da.orders@cookmedical.com Finland: +358 972519996, fi.orders@cookmedical.com France: +33 171230269, fr.orders@cookmedical.com Germany: +49 6950072804, de.orders@cookmedical.com Hungary: +36 17779199, hu.orders@cookmedical.com Ireland: +353 61239252, ie.orders@cookmedical.com Italy: +39 0269682853, it.orders@cookmedical.com Netherlands: +31 202013367, nl.orders@cookmedical.com Norway: +47 23162968, no.orders@cookmedical.com Spain: +34 912702691, es.orders@cookmedical.com Sweden: +46 858769468, se.orders@cookmedical.com Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com



USA Website: cookmedical.com EDI: cookmedical.com/edi.do Americas: Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335 E-mail: orders@cookmedical.com Australia: Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283 E-mail: cau.custserv@cookmedical.com

AI, ESC, IR, OHNS, PI, RH, SUR-A4