



An alternative for ductal access.



EchoTip[®] Ultra

HD ULTRASOUND ACCESS NEEDLE

When you have a pancreatobiliary access case that is too challenging for ERCP, you need an alternative. That's why Cook Medical created the HD Ultrasound Access Needle—the first therapeutic EUS needle indicated for gaining access. And, because it has our highly echogenic dimple pattern and ergonomic handle, it can help assure that the needle is within the target and that the procedure is stable and secure.



EchoTip® Ultra

HD ULTRASOUND ACCESS NEEDLE

This device is used to access or sample submucosal and extramural lesions of the gastrointestinal tract, access of the following: the intra- or extra-hepatic bile ducts, pancreatic ducts, cystic duct, gallbladder or for delivery of injectable materials into tissues through the accessory channel of an ultrasound endoscope.

Supplied sterile and is disposable - intended for single use only.

- Dimple pattern promotes ultrasound visibility
- Naturally contoured handle delivers precise needle control and stability
- Atraumatic needle tip designed to allow passage of a .035" wire guide to maintain access
- Beveled stylet used to puncture and gain access
- Stylet cap designed to maintain stylet position on puncture

| Order Number | Reference Part Number | Needle Gage | Needle Extension cm | Sheath Fr | Stylet Tip | Minimum Accessory Channel mm |
|--------------|-----------------------|-------------|---------------------|-----------|------------|------------------------------|
| G52012 | ECHO-HD-19-A | 19 | 0-8 | 5.2- 4.3 | Beveled | 2.8 |

5 cc and 10 cc vacuum locking syringe included.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4