


Precision* stenting
from the hilum
through the papilla.



NOW

Available with
10 cm & 12 cm
Lengths

Zilver 635[®] Biliary
SELF-EXPANDING STENT







*ZILBS-635 verification testing on file for deployment accuracy.

Warning: The safety and effectiveness of this device for use in the vascular system have not been established.

Hilar stenting with a 6 Fr introducer.

Thanks to a slim 6 Fr introducer, you can now deliver all the benefits of the Zilver to the hilar region. The Zilver 635 balances important traits - such as nonforeshortening, radial force, flexibility and patency allowing confidence in stent performance while you focus on providing better results for your patient.





Product Specifications

Order Number	Reference Part Number	Stent Diameter mm	Stent Length cm	Introducer Fr	Introducer Length cm	Wire Guide [†] Diameter inch
G50616	ZILBS-635-6-4	6	4	6	200	.035
G50619	ZILBS-635-6-6	6	6	6	200	.035
G50622	ZILBS-635-6-8	6	8	6	200	.035
G23811	ZILBS-635-6-10*	6	10	6	200	.035
G23812	ZILBS-635-6-12*	6	12	6	200	.035
G50617	ZILBS-635-8-4	8	4	6	200	.035
G50620	ZILBS-635-8-6	8	6	6	200	.035
G50623	ZILBS-635-8-8	8	8	6	200	.035
G23813	ZILBS-635-8-10*	8	10	6	200	.035
G23814	ZILBS-635-8-12*	8	12	6	200	.035
G50618	ZILBS-635-10-4	10	4	6	200	.035
G50621	ZILBS-635-10-6	10	6	6	200	.035
G50624	ZILBS-635-10-8	10	8	6	200	.035
G23815	ZILBS-635-10-10*	10	10	6	200	.035
G23816	ZILBS-635-10-12*	10	12	6	200	.035

*Check for availability. †Recommend Acrobat® 2 wire guide, sold separately. Minimum accessory channel 2.8 mm.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800772222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-A4