


Maintain access
at all times.



Fusion[®]OASIS[®]
ONE ACTION STENT INTRODUCTION SYSTEM

DUAL PLATFORM ERCP

ACCESS

STONE MANAGEMENT

STRICTURE MANAGEMENT





Fusion[®]OASIS[®]

ONE ACTION STENT INTRODUCTION SYSTEM

Used for endoscopic biliary stent placement to drain obstructed bile ducts.

- System includes a pushing catheter and guiding catheter. (The Fusion Pushing Catheter is also sold separately.)
- The integrated guiding and pushing catheters provide pushability when negotiating tortuous anatomy.
- It is compatible with a wide variety of our plastic stents.
- Eliminate the need for recannulation when placing multiple stents, repositioning a stent, or fully removing the stent from the ductal system.
- Reposition, recapture before deployment, or remove plastic stents entirely while maintaining ductal access with Fusion capability.

Order Number	Reference Part Number	Pushing Catheter Fr	Guiding Catheter Fr	Pushing/ Guiding Catheter Length cm	Wire Guide Diameter inch	Minimum Accessory Channel Diameter mm
G31526	FS-OA-8.5	8.5	5	170/195	0.035	4.2
G31527	FS-OA-10	9	5.5	170/195	0.035	4.2

Not for use with pigtail stents. 1 radiopaque band at IDE port.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu

Fusion Pushing Catheter

Order Number	Reference Part Number	Fr	Length cm	Minimum Accessory Channel Diameter mm
G31550	FS-PC-5	5	170	3.2
G31551	FS-PC-7	6	170	3.2

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4