

Assist in cell collection in the biliary system.



Fusion[®]
CYTOLOGY BRUSH

DUAL PLATFORM ERCP

ACCESS

STONE MANAGEMENT

STRICTURE MANAGEMENT

COOK[®]
MEDICAL



Fusion®

CYTOLOGY BRUSH

Used for the collection of cells in the biliary system.

- The angled bristles are designed to assist in cell collection.
- The 8 Fr cytology brush has a 1.5 cm long filiform tip.
- The Mark V System® provides endoscopic reference at 5 cm, 10 cm, and 15 cm from the distal tip of the catheter.

Order Number	Reference Part Number	Brush Length cm	Catheter Fr/Length cm	Filiform Tip Length cm	Wire Guide Diameter inch	Minimum Accessory Channel mm
G31525	FS-CB-1.5-5	2.5	8/200	1.5	0.035	4.2

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details. Please see product risk information in the IFU at cookmedical.eu

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4