


BILIARY SURGERY
DEVICE FAMILY

Expansive options for biliary stone management


A surgeon saw a need. We listened.
No matter your approach, our
combined sets and individual devices
can help you manage biliary stones.



DIAGNOSE

Access, view, and diagnose bile duct stones.

Fanelli Cholangiography Catheter Set



The integrated wire guide is designed to help you navigate the small, complex ducts.

These centimeter markings can help ensure you have accurate positioning.

This flexible tip with smooth edges is designed to help you pass the catheter into the cystic duct.

REMOVE

Our stone management devices are sold as kits and as individual pieces to fit any patient's need.




This tipless version of the stone extractor can be used with a scope.

NCompass®
NITINOL STONE EXTRACTOR



This tipped version can help you position it under fluoroscopy.

DIAGNOSE - INTRAOPERATIVE CHOLANGIOGRAPHY



Fanelli Cholangiography Catheter Set

Used for diagnostic evaluation of the bile ducts during laparoscopic cholecystectomy procedures.

The set includes a catheter with integrated wire guide, introducer sheath needle, and three-way multipurpose adapter.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Wire Guide Diameter inch	Wire Guide Length cm
G12950	C-FCC-100	5.0	65	0.021	90



Franklin Endoscopic Cholangiography Set


Used for diagnostic evaluation of the bile ducts during laparoscopic cholecystectomy procedures.

The set includes a catheter, wire guide, and stop-cock.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Wire Guide Diameter inch	Wire Guide Length cm
G10456	C-ECS-400-60-FECS	4.0	60	0.025	145

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

DIAGNOSE - INTRAOPERATIVE CHOLANGIOGRAPHY




Mixer Endoscopic Cholangiography Set

Used for diagnostic evaluation of the bile ducts during laparoscopic cholecystectomy procedures.

The set includes a catheter, needle stylet, and stop-cock.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Stylet Diameter inch	Stylet Length cm
G08045	C-OECS-400-MIXTER	4.0	50	0.025	52
G09331	C-OECS-500-MIXTER	5.0	50	0.035	52




Olsen Endoscopic Cholangiography Set

Used for access and diagnostic evaluation of the bile ducts during laparoscopic cholecystectomy procedures.

The set includes a catheter, introducer sheath needle, and stop-cock.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Introducer Sheath Needle gage	Introducer Sheath Needle Length cm
G05926	C-OECS-300	3.0	43	16	12
G05927	C-OECS-400	4.0	43	14	12



Petelin Endoscopic Cholangiography Set


Used for access and diagnostic evaluation of the bile ducts during laparoscopic cholecystectomy procedures.

The set includes a catheter, needle, and stop-cock.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Needle gage	Catheter Needle Length cm
G10664	C-OECS-400-PECS	4.0	60	14	7.5

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

REMOVE - STONE MANAGEMENT/COMMON BILE DUCT EXPLORATION



Appel-Berci Cystic Duct Introducer Set

Used for placement of choledochoscopes, cholangiography catheters, and other instruments.

The set includes an introducer sheath, dilator, needle, and wire guide.


Order Number	Reference Part Number	Sheath ID mm	Sheath Length cm	Introducer Fr	Accepts Instrumentation Fr
G08297	C-CDIS-4.0-15-BERCI	4.0	15	12.0	12.0



Roadrunner® Biliary Wire Guide

Used to assist in cannulation for laparoscopic and general surgical procedures for the dilation of the cystic duct to facilitate common bile duct exploration.

Order Number	Reference Part Number	Diameter inch	Length cm	Taper Length cm	Floppy Tip Length cm	Tip Configuration
G58996	RFSPC-35-145-CDES	0.035	145	11.5	3	Angled




Advance® Biliary Balloon Catheter

Used for the dilation of the cystic duct to facilitate common bile duct exploration.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Accepts Wire Guide Diameter inch	Inflated Balloon Diameter mm	Balloon Length cm
G58803	C-PTA5-35-80-6-4.0	5.0	80	0.035	6	4
G58804	C-PTA5-35-80-8-4.0	5.0	80	0.035	8	4

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.


REMOVE - STONE MANAGEMENT/COMMON BILE DUCT EXPLORATION



Cook Sphere® Inflation Device

Used to inflate the balloon with fluid, to monitor the pressure in the balloon during the procedure, and to deflate the balloon after dilation.


Order Number	Reference Part Number	Fluid Capacity mL	Maximum Pressure atm
G34903	CID-60-15	60	15



Three-Way Multipurpose Adapter

Used for saline flush or contrast for cholangiography.

Order Number	Reference Part Number
G13070	C-MPA-200




NCompass® Nitinol Stone Extractor

Used for extraction of stones or debris during biliary surgical procedures.

Order Number	Reference Part Number	Catheter Fr	Length cm	Basket Configuration	Deployed Basket Diameter cm
G36250	C-NTSE-2.4-115-NC3	2.4	115	3 wires, tipped	1.3
G36251	C-NTSE-2.4-115-NCT4	2.4	115	4 wires, tipless	1.3

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

REMOVE - STONE MANAGEMENT/COMMON BILE DUCT EXPLORATION




Strange Bile Duct Stone Exploration Set

Used for cholangiography and bile duct stone retrieval using fluoroscopy. The device is indicated for adults only.

The set includes a double-lumen catheter with preloaded inner catheter, wire guide, NCompass® Nitinol Stone Extractor with Tuohy-Borst adapter, and an introducer with preloaded needle.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm
G09108	C-CSRC-7.5D-40-STRANGE	7.5	40



Nathanson Cystic/Common Bile Duct Exploration Set

Used for extraction of bile duct stones via a transcystic approach.

The set includes a catheter, stone basket extractor, and Tuohy-Borst Adapter.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Basket Configuration	Basket Material	Deployed Basket Diameter cm
G03298	CDES-550-NATHANSON	5.5	70	4 wires	Stainless steel	2.0

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

No matter your approach, we have the tools to help.

DIAGNOSE

Intraoperative cholangiography

- Fanelli Cholangiography Catheter Set
- Franklin Endoscopic Cholangiography Set
- Mixer Endoscopic Cholangiography Set
- Olsen Endoscopic Cholangiography Set
- Petelin Endoscopic Cholangiography Set

REMOVE

Transcystic choledochoscopy

- Appel-Berci Cystic Duct Introducer Set
- Roadrunner® Biliary Wire Guide
- Advance® Biliary Balloon Catheter
- Cook Sphere® Inflation Device
- Three-Way Multipurpose Adapter
- NCompass® Nitinol Stone Extractor

Transcystic fluoroscopy

- Appel-Berci Cystic Duct Introducer Set
- Roadrunner® Biliary Wire Guide
- Advance® Biliary Balloon Catheter
- Cook Sphere® Inflation Device
- Three-Way Multipurpose Adapter
- NCompass® Nitinol Stone Extractor
- Strange Bile Duct Stone Exploration Set
- Nathanson Cystic/Common Bile Duct Exploration Set

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4