

Take control of your tissue acquisition and foreign-body retrieval procedures.



TISSUE MANAGEMENT

TARGET | **ACQUIRE** | TREAT

Captura Pro™

BIOPSY FORCEPS

Captura Pro Biopsy Forceps are available in various sizes and configurations to meet your clinical needs and preferences. Choose the right cup size, shaft length, and forceps type for the procedure.





Captura Pro™

BIOPSY FORCEPS

Used to obtain endoscopic mucosal tissue biopsies and/or for foreign-body retrieval.

- The dual drive-wire design enables the biopsy cups to pivot.
- The handle is designed for shorter throw length.
- The markings are designed to signal position and movement.

Order Number	Reference Part Number	Cup Diameter mm	Shaft Length cm	Forceps Type
Standard Forceps				
G47687	BF-2.2S-160-20	2.2	160	nonspiked
G47688	BF-2.2S-160-40	2.2	160	nonspiked
G47689	BF-2.2S-160SP-20	2.2	160	spiked
G47686	BF-2.2S-160SP-40	2.2	160	spiked
G47690	BF-2.2S-230-20	2.2	230	nonspiked
G47691	BF-2.2S-230-40	2.2	230	nonspiked
G47692	BF-2.2S-230SP-20	2.2	230	spiked
G47693	BF-2.2S-230SP-40	2.2	230	spiked
Large Forceps				
G47694	BF-2.4SL-160-20	2.4	160	nonspiked
G47695	BF-2.4SL-160-40	2.4	160	nonspiked
G47696	BF-2.4SL-160SP-20	2.4	160	spiked
G47697	BF-2.4SL-160SP-40	2.4	160	spiked
G47698	BF-2.4SL-230-20	2.4	230	nonspiked
G47699	BF-2.4SL-230-40	2.4	230	nonspiked
G50696	BF-2.4SL-230SP-20	2.4	230	spiked
G50697	BF-2.4SL-230SP-40	2.4	230	spiked

Order Number	Reference Part Number	Cup Diameter mm	Shaft Length cm	Forceps Type
Jumbo Forceps				
G56992	BF-2.8SJ-160-20	2.8	160	nonspiked
G56993	BF-2.8SJ-160-40	2.8	160	nonspiked
G56994	BF-2.8SJ-160SP-20	2.8	160	spiked
G56995	BF-2.8SJ-160SP-40	2.8	160	spiked
G56996	BF-2.8SJ-230-20	2.8	230	nonspiked
G56997	BF-2.8SJ-230-40	2.8	230	nonspiked
G56998	BF-2.8SJ-230SP-20	2.8	230	spiked
G56999	BF-2.8SJ-230SP-40	2.8	230	spiked
Max Forceps				
G57000	BF-3.3SM-230-20	3.3	230	nonspiked
G57001	BF-3.3SM-230-40	3.3	230	nonspiked
G57002	BF-3.3SM-230SP-20	3.3	230	spiked
G57003	BF-3.3SM-230SP-40	3.3	230	spiked
Captura® Mini Biopsy Forceps				
G47742	DBF-1.8-S	1.8	110	nonspiked
G53006	DBF-1.8-160-S	1.8	160	nonspiked
Captura® Hot Biopsy Forceps				
G31583	HDBF-2.4-230-S	2.4	230	nonspiked
G56135	HDBF-2.4SN-230-S	2.4	230	nonspiked

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Please refer to the product's Instructions for Use (IFU) for full prescribing information, warnings, precautions, contraindications, and potential adverse events.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4