

Deliver a fine staining-solution mist.



Glo-Tip[®]
SPRAY CATHETER

When spraying the gastrointestinal mucosa to identify suspicious lesions, a fine mist is essential to optimally stain the tissue. The Glo-Tip Spray Catheter delivers a fine mist, allowing you to perform effective procedures. The PTFE sheath provides low coefficient of friction for an even flow of the staining solution and for ease of insertion.



Glo-Tip®

SPRAY CATHETER

Used in conjunction with an endoscope to spray the gastrointestinal mucosa.

It is supplied sterile and is disposable. It is intended for single use only.

- The PTFE sheath provides a low coefficient of friction for ease of insertion, allowing for an even flow of staining solution.
- The radiopaque tip enhances fluoroscopic visualization, if needed.
- It has a smooth, metal spray tip.
- Its disposability eliminates cleaning and reprocessing costs and saves time.

Order Number	Reference Part Number	Catheter Fr/mm	Catheter Length cm	Stylet Length cm
G24892	GT-7-SPRAY	6.9/2.3	235	100

Minimum accessory channel diameter 2.8 mm.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com-A4



AI, ESC, IR, OHNS, PI, RH, SUR-A4