

Simple, secure wire management.



SnapLoc™

WIRE GUIDE LOCKING DEVICE

During your next ERCP procedure, consider using SnapLoc to securely lock the wire guide in place and eliminate the need for guess work. The simple, one-step locking action gives you confidence in a firmly fastened attachment while maintaining access and during exchanges.



SnapLoc™

WIRE GUIDE LOCKING DEVICE

This device is an accessory to be used with endoscopic biliary and pancreatic devices to lock the wire guide(s) in place and to prevent reflux of bodily fluids during ERCP.

- Allows you to manage up to two wires at a time.
- Helps reduce reflux of bodily fluids during device exchange.¹
- Provides a one-step locking action for confident attachment.
- Facilitates less clutter at the instrument port because of the low-profile, intuitive design.

Order Number	Reference Part Number	Endoscope Manufacturer	Box Quantity
G59683	SNAPLOC-O	Olympus®	10
G59686	SNAPLOC-P	Pentax®	10
G59685	SNAPLOC-F	Fujifilm®	10

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Refer to the Instructions for Use (IFU) for complete product information and full prescribing information including indications, contraindications, warnings, precautions, and potential adverse events. Please see product risk information in the IFU at www.cookmedical.com

OLYMPUS is a registered trademark of Olympus Corporation (Tokyo, Japan), Olympus America Inc. or their affiliated entities.

PENTAX is a registered trademark of Pentax Imaging Company.

FUJIFILM is a registered trademark of FUJIFILM corporation.

¹Internal Cook documents: #D00266104 and #D00266105.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4