


Maintain access
at all times.



Fusion[®]OASIS[®]

ONE ACTION STENT INTRODUCTION SYSTEM

DUAL PLATFORM ERCP

ACCESS

STONE MANAGEMENT

STRICTURE MANAGEMENT



Fusion[®]OASIS[®]

ONE ACTION STENT INTRODUCTION SYSTEM

Used for endoscopic biliary stent placement to drain obstructed bile ducts.

- System includes a pushing catheter and guiding catheter. (The Fusion Pushing Catheter is also sold separately.)
- The integrated guiding and pushing catheters provide pushability when negotiating tortuous anatomy.
- It is compatible with a wide variety of our plastic stents.
- Eliminate the need for recannulation when placing multiple stents, repositioning a stent, or fully removing the stent from the ductal system.
- Reposition, recapture before deployment, or remove plastic stents entirely while maintaining ductal access with Fusion capability.

Order Number	Reference Part Number	Pushing Catheter Fr	Guiding Catheter Fr	Pushing/ Guiding Catheter Length cm	Wire Guide Diameter inch	Minimum Accessory Channel Diameter mm
G31526	FS-OA-8.5	8.5	5	170/195	0.035	4.2
G31527	FS-OA-10	9	5.5	170/195	0.035	4.2

Not for use with pigtail stents. 1 radiopaque band at IDE port.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Fusion Pushing Catheter

Order Number	Reference Part Number	Fr	Length cm	Minimum Accessory Channel Diameter mm
G31550	FS-PC-5	5	170	3.2
G31551	FS-PC-7	6	170	3.2

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4