


BILIARY SURGERY
DEVICE FAMILY

Options for biliary stone management

Our combined sets and individual devices
can help you manage biliary stones.




DIAGNOSE

Access, view, and diagnose bile duct stones.

Mixer Endoscopic Cholangiography Set

- The full set allows for simple introduction, is designed for ease of use, and creates efficiencies within your surgical suite.
- The pre-loaded needle stylet has a diamond-shaped tip for ease of percutaneous introduction (without the use of added access sheath or port).
- The unique combination of memory and pushability is built into the catheter.




REMOVE


Our stone management devices are sold as individual pieces to fit any patient's need.

NCompass[®]
NITINOL STONE EXTRACTOR

This tipless version of the stone extractor can be used with a scope.

This tipped version can help you position the extractor under fluoroscopy.






Mixer Intraoperative Cholangiography Set

Used for diagnostic evaluation of the bile ducts during laparoscopic cholecystectomy procedures.

The set includes a catheter, needle stylet, and stopcock.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Stylet Diameter inch	Stylet Length cm
G08045	C-OECS-400-MIXTER	4.0	50	0.025	52
G08045	C-OECS-500-MIXTER	5.0	50	0.035	52

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.




Olsen Intraoperative Cholangiography Set

Used for access and diagnostic evaluation of the bile ducts during laparoscopic cholecystectomy procedures.

The set includes a catheter, introducer sheath needle, and stopcock.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Introducer Sheath Needle gage	Introducer Sheath Needle Length cm
G05926	C-OECS-300	3.0	43	16	12
G05926	C-OECS-400	4.0	43	14	12

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.



Cystic Duct Introducer Set

Used for placement of choledochoscopes, cholangiography catheters, and other instruments.

The set includes an introducer sheath, dilator, needle, and wire guide.

Order Number	Reference Part Number	Sheath ID mm	Sheath Length cm	Introducer Fr	Accepts Instrumentation Fr
G08297	C-CDIS-4.0-15-BERCI	4.0	15	12.0	12.0

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.




Roadrunner® Biliary Wire Guide

Used to assist in cannulation for laparoscopic and general surgical procedures for the dilation of the cystic duct to facilitate common bile duct exploration.

Order Number	Reference Part Number	Diameter inch	Length cm	Taper Length cm	Floppy Tip Length cm	Tip Configuration
G58996	RFSPC-35-145-CDES	0.035	145	11.5	3	Angled

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.



Advance® Biliary Balloon Catheter

Used for the dilation of the cystic duct to facilitate common bile duct exploration.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Accepts Wire Guide Diameter inch	Inflated Balloon Diameter mm	Balloon Length cm
G58803	C-PTA5-35-80-6-4.0	5.0	80	0.035	6	4
G58804	C-PTA5-35-80-8-4.0	5.0	80	0.035	8	4

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.



Cook Sphere® Inflation Device

Used to inflate the balloon with fluid, to monitor the pressure in the balloon during the procedure, and to deflate the balloon after dilation.

Order Number	Reference Part Number	Fluid Capacity mL	Maximum Pressure atm
G34903	CID-60-15	60	15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.




Multiport Adapter for Cholangiography

Used for saline flush or contrast for cholangiography.

Order Number	Reference Part Number
G13070	C-MPA-200

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.




NCompass® Nitinol Stone Extractor

Used for extraction of stones or debris during biliary surgical procedures.

Order Number	Reference Part Number	Catheter Fr	Length cm	Basket Configuration	Deployed Basket Diameter cm
G36250	C-NTSE-2.4-115-NC3	2.4	115	3 wires, tipped	1.3
G36251	C-NTSE-2.4-115-NCT4	2.4	115	4 wires, tipped	1.3

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.



Nathanson Cystic/Common Bile Duct Exploration Set

Used for extraction of bile duct stones via a transcystic approach.

The set includes a catheter, stone basket extractor, and Tuohy-Borst Adapter.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Basket Configuration	Basket Material	Deployed Basket Diameter cm
G03298	CDES-550-NATHANSON	5.5	70	4 wires	Stainless steel	2.0

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.


No matter your approach, we have the tools to help.

<p>DIAGNOSE</p>	<p>Intraoperative cholangiography</p> <ul style="list-style-type: none"> • Mixer Intraoperative Cholangiography Set • Olsen Intraoperative Cholangiography Set
<p>REMOVE</p>	<p>Transcystic choledochoscopy</p> <ul style="list-style-type: none"> • Cystic Duct Introducer Set • Roadrunner® Biliary Wire Guide • Advance® Biliary Balloon Catheter • Cook Sphere® Inflation Device • Multiport Adapter for Cholangiography • NCompass® Nitinol Stone Extractor
	<p>Transcystic fluoroscopy</p> <ul style="list-style-type: none"> • Cystic Duct Introducer Set • Roadrunner® Biliary Wire Guide • Advance® Biliary Balloon Catheter • Cook Sphere® Inflation Device • Multiport Adapter for Cholangiography • NCompass® Nitinol Stone Extractor • Nathanson Cystic/Common Bile Duct Exploration Set

Customer Service

EU website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, emea.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 21612723770, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com
 EDI: cookmedical.com/edi.do
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 Email: customersupport@cookmedical.com
Australia:
 Phone: +61 734346000, 1800772222, Fax: +61 734346001, 1800077283
 Email: cau.custserv@cookmedical.com



AI, ESC, IR, PI-8.5X11