

Helping to meet your specific clinical needs



Tri-Ex[®]

MULTIPLE SIZE EXTRACTION BALLOON

The Tri-Ex Multiple Size Extraction Balloon delivers inflation diameters of 8.5, 12, or 15 mm to meet your clinical needs. And, because there are three lumens, you can place a 0.035 inch wire guide and inject contrast during stone extraction.



Tri-Ex®

MULTIPLE SIZE EXTRACTION BALLOON

Used for endoscopic removal of biliary stones.

- Reference marks on the syringe allow for three inflation diameters (8.5, 12, or 15 mm) to meet your specific clinical needs
- Three separate lumens allow for contrast injection without compromising wire guide movement
- Contrast injection sites available above or below the balloon
- Radiopaque catheter enhances fluoroscopic visualisation
- Radiopaque bands identify balloon location fluoroscopically
- Black etched mark provides endoscopic identification that the balloon is about to exit the papilla

Order Number	Reference Part Number	Catheter Size Fr	Catheter Length cm	Balloon Diameter mm	Wire Guide Diameter inch	Injection Port	Minimum Accessory Channel mm
G26856	TXR-8.5-12-15-A	6.9–5.4	200	8.5 12 15	0.035	above balloon	2.8
G26857	TXR-8.5-12-15-B	6.9–5.4	200	8.5 12 15	0.035	below balloon	2.8
G50695	TXR-ME	6.6	275	8.5 12 15	0.035	above balloon	3.2

Caution: This product contains natural rubber latex, which may cause allergic reactions.

Some products or part numbers may not be available in all markets.

Contact your local Cook Medical representative or Customer Support & Distribution for details.

Refer to the Instructions for Use (IFU) for complete product information and full prescribing information including indications, contraindications, warnings, precautions, and potential adverse events.

Please see product risk information in the IFU at www.cookmedical.eu.

This device is supplied sterile and is disposable. Intended for single use only.

Customer Service

EU Website: cookmedical.eu

EDI: cookmedical.eu/edi

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 972519996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com-A4



AI, ESC, IR, OHNS, PI, RH, SUR-A4