

# Streamline your stone management procedures.



## Fusion Quattro<sup>®</sup>

EXTRACTION BALLOON

DUAL PLATFORM ERCP

ACCESS

STONE MANAGEMENT

STRICTURE MANAGEMENT



## Fusion Quattro®

### EXTRACTION BALLOON

Used for endoscopic removal of stones in the biliary system and for contrast injection.

- 4 inflation diameters in 1 device to meet your specific clinical needs.
- Ease of injection with contrast exiting the distal tip for above-balloon injections.
- Radiopaque bands enhance fluoroscopic visualisation of the extraction balloon.
- Reference marks on syringe indicate each inflation diameter.

Order Number	Reference Part Number	Balloon Diameter mm	Catheter Size Fr	Catheter Length cm	Injection Port	Wire Guide Diameter inch
G48225	FS-QEB-A	8.5, 10, 12, 15	6.6	200	Above Balloon	0.035
G48226	FS-QEB-B	8.5, 10, 12, 15	6.6	200	Below Balloon	0.035
G31921	FS-QEB-XL-A	12, 15, 18, 20	6.6	200	Above Balloon	0.035
G31922	FS-QEB-XL-B	12, 15, 18, 20	6.6	200	Below Balloon	0.035

Minimum accessory channel 3.2 mm.

**Caution:** This product contains natural rubber latex, which may cause an allergic reaction.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Refer to the Instructions for Use (IFU) for complete product information and full prescribing information including indications, contraindications, warnings, precautions, and potential adverse events.

Please see product risk information in the IFU at [www.cookmedical.eu](http://www.cookmedical.eu)

### Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Iceland: +354 800 7615, [is.orders@cookmedical.com](mailto:is.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)

EDI: [cookmedical.com/edi](http://cookmedical.com/edi)

#### Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [customersupport@cookmedical.com](mailto:customersupport@cookmedical.com)

#### Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



AI-ESC-IR-OHNS-PI-RH-SUR-A4