

Reliability matters.

Target and obtain tissue with a needle designed to deliver consistency across each pass.



EchoTip® ClearCore™

EUS BIOPSY NEEDLE

Smooth puncture

The needle offers 45% less puncture force and shows consistency across individual passes compared to the leading 19 gage FNB needle. The Franseen tip is made from a cobalt-chromium (CoCr) alloy for durability and to minimize deformation.

Enhanced control

The ergonomic, contoured handle facilitates needle control and stability. The nitinol stylet supports needed puncture force during targeting and sampling. And the echogenic dimple pattern on the needle promotes visibility during the procedure.

Targeted precision

The PEEK (polyetheretherketone) sheath promotes flexibility and ease of use when targeting lesions in tortuous anatomy, limiting the need for scope adjustments even in difficult positions.



EchoTip® ClearCore™

EUS BIOPSY NEEDLE

The **EchoTip ClearCore EUS Biopsy Needle** is used with an ultrasound endoscope for fine needle biopsy (FNB) of submucosal and extramural lesions, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the gastrointestinal tract.

The needle is compatible with a range of endoscopes.

It is supplied sterile and disposable. It is intended for single use only.

Order Number	Reference Part Number	Sheath Fr	Needle gage	Needle Tip Bevel	Adjustable Needle Extension length cm	Minimum Accessory Channel Diameter mm
G59742	ECHO-BX-19	5.4	19	Franseen	0-8	2.8

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-8.5X11