


Flexibility matters.

Target and acquire tissue with a needle designed to maximise flexibility and minimise puncture force.



EchoTip® ClearCore™ EUS BIOPSY NEEDLE

Enhanced flexibility

The spring-coiled sheath—coated inside and out—is designed to maximise flexibility for accessing difficult anatomy^{1,2} and gives you exceptional control to advance the needle without adjusting your scope.³

Optimised puncture

The echogenic Franseen needle—made from a cobalt-chromium (CoCr) alloy to minimise deformation⁴—allows you to puncture tissue with minimum force^{3,5} to rapidly target lesions.

Controlled precision

The ergonomic, contoured handle facilitates needle control and stability, and the ReCoil® stylet's unique design allows you to easily manage the stylet during the procedure.



EchoTip® ClearCore™

EUS BIOPSY NEEDLE

The **EchoTip ClearCore EUS Biopsy Needle** is used with an ultrasound endoscope for fine needle biopsy (FNB) of submucosal and extramural lesions, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the gastrointestinal tract.

The needle is compatible with a range of endoscopes.²

It is supplied sterile and disposable. It is intended for single use only.

Order Number	Reference Part Number	Needle gauge	Sheath Fr	Needle Extension Length cm	Needle Tip Bevel	Minimum Accessory Channel Diameter mm
G60992	ECHO-BX-22	22	5.3	0-8	Franseen	2.8

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.

Data on file:

1. D00527893, D00516367, D00516371, D00521988
2. D00527893, D00516367, D00516371
3. D00537459
4. Keehan E, Gergely L. Catheter and specialty needle alloys. Poster session presented at: Materials & Processes for Medical Devices Conference & Exposition; August 2009; Minneapolis, MN.
5. D00300922, D00527893, D00516367, D00516371, D00525720

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
 Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com
 Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-8.5X11