

Sizing up your polypectomy procedure.



AcuSnare[®]
POLYPECTOMY SNARE

Polyps come in a wide variety of shapes and sizes. So do our AcuSnares. Available in multiple configurations—and ranging in flexibility—the AcuSnares allow you to make the best choice to meet the procedural challenges you face. Cook Medical AcuSnares. Good choice.





AcuSnare®

POLYPECTOMY SNARE

This device is used endoscopically in the removal and cauterization of sessile polyps and pedunculated polyps from within the gastrointestinal tract.

- Braided stainless steel snare for smooth and concise excision
- 6.9 FR, 240 cm PTFE sheath delivers low coefficient of friction for smooth insertion

Order Number	Reference Part Number	Snare Size cm	Sheath French Size Fr	Sheath Length cm	Snare Configuration	Box Quantity
AcuSnare®						
G22631	ASM-1-S	1.5 x 3	6.9	240	Mini Oval	10
G22629	AS-1-S	2.5 x 5.5	6.9	240	Standard Oval	10
G22630	ASJ-1-S	3 x 6	6.9	240	Jumbo Oval	10
G22701	ASMH-1-S	1.5 x 2.5	6.9	240	Mini Hexagonal	10
G22700	ASH-1-S	3 x 4.5	6.9	240	Hexagonal	10
G22647	ASDB-15-015-S	1.5	6.9	240	Duck-Bill	10
G22649	ASDB-25-015-S	2.5	6.9	240	Duck-Bill	10
Soft AcuSnare®						
G22704	SASMM-1-S	1 x 1.5	6.9	240	Micro Mini Oval	10
G22633	SASM-1-S	1.5 x 3	6.9	240	Mini Oval	10
G22632	SAS-1-S	2.5 x 5.5	6.9	240	Standard Oval	10
G22703	SASMH-1-S	1.5 x 2.5	6.9	240	Mini Hexagonal	10
G22702	SASH-1-S	3 x 4.5	6.9	240	Hexagonal	10

Minimum accessory channel 2.8 mm. Active cord available separately.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.

Customer Service

EU website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, emea.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 21612723770, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

Email: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

Email: cau.custserv@cookmedical.com



AI, ESC, IR, PI-A4