


# Directly measure portal pressure under EUS.



**EchoTip Insight<sup>®</sup>**  
PORTOSYSTEMIC PRESSURE GRADIENT




## EchoTip Insight®

PORTOSYSTEMIC PRESSURE GRADIENT

**EchoTip Insight Needle:** This device is used to access the hepatic and portal veins or their branches, via a transgastric-transhepatic approach or transduodenal-transhepatic approach at the duodenal bulb, through the accessory channel of an ultrasound endoscope to allow for direct pressure measurement.

**Connecting Tube:** The tube is used for the transfer of liquids between the EchoTip Insight components.



### Compass UniversalHG


- Digital display of the pressure reading
- Self-calibrating transducer
- Compact size

### Connecting tube

- 90 cm long connecting tube

### Endoscopic ultrasound needle

- Echogenic tip – to ensure that the needle tip is within the target
- Small gauge needle (25 gauge) to facilitate access of the portal and hepatic veins



Order Number	Reference Part Number	Needle gauge	Sheath Fr	Adjustable Needle Extension Length cm	Connecting Tube Length cm	Pressure Range mm Hg	Minimum Accessory Channel Diameter mm
G58000	ECHO-PPG-I	25	5.2	8	90	-199 to +999	2.8

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at [cookmedical.eu](http://cookmedical.eu).

### Customer Service

EU website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [emea.distributors@cookmedical.com](mailto:emea.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 21612723770, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Iceland: +354 8007615, [is.orders@cookmedical.com](mailto:is.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA website: [cookmedical.com](http://cookmedical.com)  
 EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)

**Americas:**  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 Email: [customersupport@cookmedical.com](mailto:customersupport@cookmedical.com)

**Australia:**  
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283  
 Email: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



AI, ESC, IR, PI-A4