

# A different approach to hemostasis.

For nonvariceal upper and lower GI bleeds



Initial bleed



Hemospray application



Follow-up

**Hemospray<sup>®</sup>**  
ENDOSCOPIC HEMOSTAT

For more information on Hemospray,  
visit [hemospray.cookmedical.com](http://hemospray.cookmedical.com).



Images are courtesy of Dr. John Morris,  
Glasgow Royal Infirmary, Glasgow, Scotland.



## Hemospray®

ENDOSCOPIC HEMOSTAT

This device is used for hemostasis of nonvariceal gastrointestinal bleeding.

This product is supplied sterile and is disposable. It is intended for single use only.

- Nonspecific targeting: Powder is sprayed toward the source of the bleed, so it does not require the en face approach and the precise targeting of other modalities.
- Nonthermal: No tissue changes occur as sometimes experienced with thermal modalities.
- Nontraumatic: Since no direct force is applied at the treatment site, the powder minimizes the risk of tissue trauma as experienced with other modalities.
- Noncontact: Aerosol delivery system eliminates the need for direct mechanical or contact forces with the treatment site by the delivery system.
- CO<sub>2</sub>-powered delivery system: CO<sub>2</sub> helps minimize postprocedural pain compared to air.<sup>1</sup>

Order Number	Reference Part Number	Catheter Length cm	Catheter Fr	Minimum Accessory Channel Diameter mm
G56572	HEMO-7	220	7.2	2.8
G21049	HEMO-10	220	9.6	3.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Please refer to the product's Instructions for Use (IFU) for full prescribing information, warnings, precautions, contraindications, and potential adverse events.

1. ASGE Technology Committee; Lo SK, Fujii-Lau LL, Enestvedt BK, et al. The use of carbon dioxide in gastrointestinal endoscopy. *Gastrointest Endosc.* 2016;83(5):857-865.

### Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Iceland: +354 800 7615, [is.orders@cookmedical.com](mailto:is.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)

EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)

#### Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

#### Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



AI-ESC-IR-OHNS-PI-RH-SUR-A4