Simple, secure wire management.



SnapLoc* WIRE GUIDE LOCKING DEVICE

During your next endoscopic retrograde cholangiopancreatography (ERCP) procedure, consider using SnapLoc to securely lock the wire guide in place and eliminate the need for guesswork. The simple, one-step locking action gives you confidence in a firmly fastened attachment while maintaining access and during exchanges.





SnapLoc™ WIRE GUIDE LOCKING DEVICE

This device is an accessory to be used with endoscopic biliary and pancreatic devices to lock the wire guide(s) in place and to prevent reflux of bodily fluids during ERCP.

- SnapLoc allows you to manage up to two wires at a time.
- It helps reduce the reflux of bodily fluids during device exchange.¹
- It provides a one-step locking action for confident attachment.
- It also helps minimise clutter at the instrument port because of the low-profile, intuitive design.

Order Number	Reference Part Number	Endoscope Manufacturer	Box Quantity
G59683	SNAPLOC-O	Olympus®	10
G59686	SNAPLOC-P	Pentax®	10
G59685	SNAPLOC-F	Fujifilm®	10

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Refer to the Instructions for Use (IFU) for complete product information and full prescribing information including indications, contraindications, warnings, precautions, and potential adverse events. Please see product risk information in the IFU at cookmedical.eu.

OLYMPUS is a registered trademark of Olympus Corporation (Tokyo, Japan), Olympus America Inc. or their affiliated entities.

PENTAX is a registered trademark of Pentax Imaging Company.

FUJIFILM is a registered trademark of FUJIFILM corporation.

1. Internal Cook documents: #D00266104 and #D00266105.

Customer Service

EU Website: cookmedical.eu EDI: cookmedical.eu/edi

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com Belgium: +32 27001702, be.orders@cookmedical.com Denmark: +45 38487607, da.orders@cookmedical.com Finland: +358 972519996, fi.orders@cookmedical.com France: +33 171230269, fr.orders@cookmedical.com Germany: +49 6950072804, de.orders@cookmedical.com Hungary: +36 17779199, hu.orders@cookmedical.com Iceland: +354 800 7615, is.orders@cookmedical.com Ireland: +353 61239252, ie.orders@cookmedical.com Italv: +39 0269682853, it.orders@cookmedical.com Netherlands: +31 202013367, nl.orders@cookmedical.com Norway: +47 23162968, no.orders@cookmedical.com Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335 E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4