


Get closure.

One-to-one rotation for precise placement and a tissue-grasping design for secure closure.



Instinct Plus[®]
ENDOSCOPIC CLIPPING DEVICE





Instinct Plus[®]

ENDOSCOPIC CLIPPING DEVICE

The Instinct Plus Endoscopic Clipping Device is endoscopically placed within the gastrointestinal tract to apply mechanical pressure to tissue and is used for endoscopic clip placement within the gastrointestinal tract for the purpose of:

- Endoscopic marking
- Haemostasis for:
 - Mucosal/submucosal defects less than 3 cm
 - Polyps less than 1.5 cm in diameter
 - Bleeding ulcers
 - Diverticula in the colon
 - Arteries less than 2 mm
 - Prophylactic clipping to reduce the risk of delayed bleeding post lesion resection
- Anchoring to affix jejunal feeding tubes to the wall of the small bowel
- As a supplementary method for closure of GI tract luminal perforations less than 20 mm that can be treated conservatively
- Anchoring to affix fully covered oesophageal self-expanding metal stents to the wall of the oesophagus in patients with fistulas, leaks, perforations, or disunion

Order Number	Reference Part Number	Clip Opening Width mm	Catheter Fr	Catheter Length cm	Quantity	MR Status	Minimum Forward Viewing Channel Diameter mm	Minimum Side Viewing Channel Diameter mm
G58010	INSC-P-7-230-S	16	7	230	10	MR Conditional	2.8	4.2

Please follow any step-by-step instructions carefully, including the risk, technical information, Device Name and intended use provided in the IFU at cookmedical.eu. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

Customer Service

EU website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, emea.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 2161 272 3770, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com
 EDI: cookmedical.com/edi.do

Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 Email: customersupport@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 Email: cau.custserv@cookmedical.com



AI, ESC, IR, PI-A4