

From cytology to histology:
Continuing to increase
your predictive values.




EchoTip ProCore[®]
HD ULTRASOUND BIOPSY NEEDLE




Get to the core issue throughout the entire GI tract.

EchoTip ProCore® allows for collection of histology with cytology while providing access to the entire GI tract, where you can target smaller lesions and increase your yields, all while decreasing your needle passes.¹




1. Kamata K, Kitano M, Yasukawa S, et al. Histologic diagnosis of pancreatic masses using 25-gauge endoscopic ultrasound needles with and without a core trap: a multicenter randomized trial. *Endoscopy*. 2016;48(7):632-638.




ReCoil® stylet helps you easily manage the stylet, minimizing the risk for contamination

Naturally contoured handle delivers precise needle control and stability


More yield, fewer needle passes.



Cell Block



Core



Coiled sheath facilitates needle flexibility (ECHO-HD-3-20-C)


Ultrasound image courtesy of Kenneth Chang, MD
H. H. Chao Comprehensive Digestive Disease Center,
University of California-Irvine.

Pathology images courtesy of James Farrell, MD
Interventional Endoscopy and Pancreatic Diseases Section
of Digestive Diseases, Yale University School of Medicine.



"Cytology is not enough in the diagnosis and treatment of all lesions. You need sufficient quality of material to correctly characterise tumours and conditions."

**Marc Giovannini, MD, PhD, Endoscopic Unit,
Paoli-Calmettes Institute Marseille, France**



Safety lock ring easily slides and locks at the desired extension

Integrated sheath adjuster offers compatibility with a full range of echoendoscopes

Metal Luer lock for secure anchoring to the echoendoscope

Product image features EchoTip ProCore 20 gauge.


Leading innovation in EUS since 1998.

Over the years, you've told us what matters to you, and we've listened. Incorporating extensive feedback provided by clinicians, we offer you our latest developments in needle flexibility and sampling capability, tailored to you and your pathologist's needs and heavily supported by published clinical data.

19 gauge needle



20 gauge needle



25 gauge needle



Product Specifications

This device is used with an ultrasound endoscope for fine needle biopsy (FNB) of submucosal lesions or tissue, mediastinal masses, lymph nodes, and intraperitoneal masses within or adjacent to the gastrointestinal tract.

Supplied sterile and is disposable – intended for single use only.

Order Number	Reference Part Number	Needle gauge	Sheath Fr	Adjustable Needle Extension cm	Stylet Tip	Minimum Accessory Channel mm
G53585	ECHO-HD-19-C	19	4.8	0-8	recessed ball tip	2.8
G34785	ECHO-HD-3-20-C*	20	7.0-7.95	0-8	recessed ball tip	3.7
G55738	ECHO-HD-25-C	25	5.2	0-8	beveled	2.8

*ReCoil stylet is only available on EchoTip ProCore 20 gauge.

Please follow any step-by-step instructions carefully, including the risk, technical information, Device Name and intended use provided in the IFU at cookmedical.eu. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

Customer Service

EU website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, emea.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 2161 272 3770, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 8007615, is.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

Email: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

Email: cau.custserv@cookmedical.com



AI, ESC, IR, PI-A4