



Chait

PERCUTANEOUS CECOSTOMY CATHETER




Give bowel control back to your patients.



Chait

PERCUTANEOUS CECOSTOMY CATHETER



- Placing the catheter is a minimally invasive, reversible procedure.
- The catheter facilitates antegrade colonic cleansing for patients who have fecal incontinence.
- The unique, low-profile trapdoor fitting is designed to lie flat against the skin.
- The catheter enables patients to manage their fecal incontinence themselves.
- The catheter is 10.2 Fr, features three sideports, and is available in three lengths.

DAWSON-MUELLER MULTIPURPOSE DRAINAGE CATHETER

The drainage catheter is placed during initial cecostomy as a temporary catheter, allowing the tract to mature. After approximately six weeks, this catheter is exchanged for the Chait Percutaneous Cecostomy Catheter.



- The Mac-Loc® locking loop helps prevent catheter displacement.
- Soft, smooth Ultrathane® material helps maintain patient comfort.
- The hydrophilic coating eases insertion.
- The catheter is 8.5 Fr and 25 cm long and is packaged with a trocar stylet, a flexible stiffener, and an introduction stiffening cannula with loading obturator.

CHAIT ENTEROSTOMY SUTURE ANCHOR SET

Components

- Introducer needle, 18 gage, 7 cm long
- Loading cartridge containing two suture anchors
- Amplatz Extra-Stiff Wire Guide, .035 inch (.089 mm) in diameter, PTFE coated, 80 cm long
- Coons dilator, 8 Fr, radiopaque, 20 cm long



Order Number	Reference Part Number	Catheter Fr	Tract Length cm
Chait Percutaneous Cecostomy Catheter			
G11632	TDCS-100	10.2	0-6
G36392	TDCS-100-M	10.2	3-9
G12685	TDCS-100-L	10.2	6-14
Chait Access Adapter with Connecting Tube*			
G11630	CTU10.0-18-CAA	10.2	-
Dawson-Mueller Multipurpose Drainage Catheter			
G09705	ULT8.5-38-25-P-5S-CLDM-HC	8.5	-
Chait Enterostomy Suture Anchor Set			
G12230	GIAS-100-CHAIT	-	-
Amplatz Extra-Stiff Wire Guide			
G04003	THSF-35-80-AES	-	-
Coons Dilator			
G03927	JCD8.0-38-20-COONS	8.0	-

*Packaged in sets of 5.

Customer Service

EMEA: EDI - cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

cookmedical.eu

Americas: EDI - cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com

