


Doctor choice, patient comfort



PERCUTANEOUS ABSCESS AND MULTIPURPOSE DRAINAGE PRODUCTS

Cook Medical offers one of the widest assortments of drainage products available, while focusing relentlessly on patient comfort.

In this catalogue you will find:

- the right choice of multipurpose drainage catheters to suit each procedure,
- the technology that only a leader can provide, including Ultrathane® construction for patient comfort, hydrophilic coating for ease of insertion, radiopacity for precise placement,
- and the product qualities that you have come to expect, such as kink-resistance and consistent hub sizes.

This catalogue focuses on our percutaneous abscess and multipurpose drainage products.

To see our full range of percutaneous drainage products, please request these additional catalogues from your Cook representative:

- **Percutaneous Biliary Intervention Products**
- **Percutaneous Urinary Drainage Products**
- **Percutaneous Access, Targeting and Accessory Products**

Be sure to have the right drainage devices for each procedure, including your most challenging and complex cases.

The choice is yours.

Contents

Non-locking Multipurpose Drainage

Multipurpose Pigtail Drainage Catheter Needle Set	4
Universal Curved Drainage Catheter	5
Non-Locking Multipurpose Drainage Catheter	6

Locking Multipurpose Drainage


Dawson-Mueller Multipurpose Drainage Catheter	7
Multipurpose Drainage Set – Micropuncture	8
Multipurpose Drainage Set – Seldinger Technique	9
Multipurpose Drainage Catheter	10

Large Lumen Drainage

Cook-Cope Locking Loop Multipurpose Drainage Catheter	11
Thal-Quick Abscess Drainage Set	12
Gordon Large-Bore Drainage Catheter	13

Double Lumen Drainage

Ring-McLean Sump Drainage Set – Straight	14
Ring-McLean Sump Drainage Set – Curved	15
Cook-Cope Ring-McLean Sump Drainage Set – Locking Pigtail	16






Multipurpose Pigtail Drainage Catheter Needle Set

Used for percutaneous external drainage in multiple drainage applications (e.g., nephrostomy, abscess drainage, and other abdominal cavity drainage procedures).

- The catheter features a non-locking pigtail for retention.
- The trocar stylet fits into the catheter and the stiffening cannula. This assembly facilitates access by direct puncture.
- The stiffening cannula protrudes from the catheter tip in order to ease the advancement of the catheter over the wire guide.
- The catheter is supplied with a pigtail straightener.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gage
G56888	DUAN-5.0-19.5-15.0-6S-PIG	5	.038	15	6	19.5
G56889	DUAN-5.0-19.5-25.0-6S-PIG	5	.038	25	6	19.5
G56890	DUAN-6.0-19.5-25.0-6S-PIG	6	.038	25	6	19.5
G56891	DUAN-7.0-19.5-25.0-8S-PIG	7	.038	25	8	19.5
G56892	DUAN-8.3-18-25.0-8S-PIG	8.3	.038	25	8	18

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Universal Curved Drainage Catheter

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess), either by direct puncture or Seldinger access technique.

- The curved tip configuration facilitates drainage in small cavities.
- Ultrathane® material enhances patient comfort.
- The Intro-Tip™ design helps to ease catheter insertion.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during insertion.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gage
G10956	ULT6.0-35-25-P-3S-UCD-HC	6	.035	25	3	19 UTW*
G10957	ULT7.0-38-25-P-3S-UCD-HC	7	.038	25	3	18
G10958	ULT8.5-38-25-P-3S-UCD-HC	8.5	.038	25	3	18
G10959	ULT10.2-38-25-P-3S-UCD-HC	10.2	.038	25	3	18

*Ultrathin wall

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Non-Locking Multipurpose Drainage Catheter

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess), either by direct puncture or Seldinger access technique.

- The catheter features a non-locking pigtail for retention.
- Ultrathane® material enhances patient comfort.
- The catheter can be inserted by direct puncture with the trocar stylet that is included or by Seldinger access technique.
- The 8.5 Fr catheter features Intro-Tip™ design that helps to ease catheter insertion.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gage
G06483	ULT8.5-38-25-P-6S-MDC*	8.5	.038	25	6	18
G06484	ULT10.2-38-25-P-6S-MDC	10.2	.038	25	6	18
G06485	ULT12.0-38-25-P-6S-MDC	12	.038	25	6	18
G06482	ULT14.0-38-25-P-6S-MDC	14	.038	25	6	18

*Features Intro-Tip

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Dawson-Mueller Multipurpose Drainage Catheter

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess), either by direct puncture or Seldinger access technique.

- The 10 mm pigtail facilitates drainage in small cavities.
- Ultrathane® material enhances patient comfort.
- The Intro-Tip™ design helps to ease catheter insertion.
- The easy-to-use Mac-Loc® locking mechanism firmly secures the pigtail.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during catheter insertion.
- The catheter also comes in an extended length version with a 24 cm long sideport segment with a radiopaque band marking the location of the most proximal sideport.

Components




Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gage
Catheter Standard Length						
G11018	ULT5.0-18-25-P-5S-CLDM-HC	5	.018	25	5	23 UTW*
G51594	ULT6.3-35-15-P-5S-CLDM-HC	6.3	.035	15	5	19 UTW*
G51595	ULT6.3-35-25-P-5S-CLDM-HC	6.3	.035	25	5	19 UTW*
G11020	ULT7.0-35-25-P-5S-CLDM-HC	7	.035	25	5	19 UTW*
G11386	ULT8.5-38-15-P-5S-CLDM-HC	8.5	.038	15	5	18
G09705	ULT8.5-38-25-P-5S-CLDM-HC	8.5	.038	25	5	18
G11387	ULT10.2-38-15-P-5S-CLDM-HC	10.2	.038	15	5	18
G09706	ULT10.2-38-25-P-5S-CLDM-HC	10.2	.038	25	5	18
G11720	ULT12.0-38-15-P-5S-CLDM-HC	12	.038	15	5	18
G10435	ULT12.0-38-25-P-5S-CLDM-HC	12	.038	25	5	18
G10436	ULT14.0-38-25-P-5S-CLDM-HC	14	.038	25	5	18
Catheter with Extended Length Sideport Segment						
G28938	ULT6.0-35-50-P-CLDMS-24-RH**	6	.035	50		

*Ultrathin wall

**Nonstock. Only flexible stiffening cannula included.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Multipurpose Drainage Set – Micropuncture

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess).

- Catheter placement is performed using micropuncture with the Neff Percutaneous Access Set that is conveniently included.
- Ultrathane® material enhances patient comfort.
- The easy-to-use Mac-Loc® locking mechanism firmly secures the 25 mm pigtail.
- The radiopaque band identifies the proximal area of the locking loop, which contains six large oval sideports to facilitate drainage.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during catheter insertion.
- The 8.5 Fr catheter features Intro-Tip™ design that helps to ease catheter insertion.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Dilators Fr
G10049	CLM-8.5-RH-NPAS-NT*	8.5	.038	25	6	9
G10050	CLM-10.2-RH-NPAS-NT	10.2	.038	25	6	8, 11
G10261	CLM-12.0-RH-NPAS-NT	12	.038	25	6	10, 14

*Features Intro-Tip

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Multipurpose Drainage Set – Seldinger Technique

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess) by Seldinger access technique.

- The set includes everything needed for an efficient catheter placement procedure using the Seldinger technique.
- The EchoTip® needle cannula enhances visibility under ultrasound.
- Ultrathane® material enhances patient comfort.
- The easy-to-use Mac-Loc® locking mechanism firmly secures the 25 mm pigtail.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during catheter insertion.
- The radiopaque band identifies the proximal area of the locking loop, which contains six large oval sideports to facilitate drainage.
- The 8.5 Fr catheter features Intro-Tip™ design that helps to ease catheter insertion.

Components



Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Dilators Fr	Needle gage/Length cm
G24874	SMS-ULT8.5-38-25-P-6S-CLM-RH*	8.5	25	6	6, 8	18/20
G24875	SMS-ULT10.2-38-25-P-6S-CLM-RH	10.2	25	6	8, 10	18/20
G24876	SMS-ULT12.0-38-25-P-6S-CLM-RH	12	25	6	8, 10, 12	18/20
G24877	SMS-ULT14.0-38-25-P-6S-CLM-RH	14	25	6	8, 10, 12, 14	18/20

*Features Intro-Tip

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Multipurpose Drainage Catheter

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess), either by direct puncture or Seldinger access technique.

- Ultrathane® material enhances patient comfort.
- The easy-to-use Mac-Loc® locking mechanism firmly secures the 25 mm pigtail.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during catheter insertion.
- The radiopaque band identifies the proximal area of the locking loop, which contains six large oval sideports to facilitate drainage.
- The 8.5 Fr catheter features Intro-Tip™ design that helps to ease catheter insertion.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gage
G11105	ULT8.5-38-15-P-6S-CLM-RH*	8.5	.038	15	6	18
G09501	ULT8.5-38-25-P-6S-CLM-RH*	8.5	.038	25	6	18
G09765	ULT8.5-38-45-P-6S-CLM-RH*	8.5	.038	45	6	18
G11210	ULT8.5-38-60-P-6S-CLM-RH* **	8.5	.038	60	6	18
G09502	ULT10.2-38-25-P-6S-CLM-RH	10.2	.038	25	6	18
G09766	ULT10.2-38-45-P-6S-CLM-RH	10.2	.038	45	6	18
G28916	ULT10.2-38-60-P-6S-CLM-RH**	10.2	.038	60	6	18
G09503	ULT12.0-38-25-P-6S-CLM-RH	12	.038	25	6	18
G09767	ULT12.0-38-45-P-6S-CLM-RH	12	.038	45	6	18
G09504	ULT14.0-38-25-P-6S-CLM-RH	14	.038	25	6	18

*Features Intro-Tip

**Nonstock

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Cook-Cope Locking Loop Multipurpose Drainage Catheter

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess), either by direct puncture or Seldinger access technique.

- The large-bore catheter features a large lumen for higher flow rate and facilitates drainage of thick fluids that will not adequately drain through smaller devices.
- Ultrathane® material enhances patient comfort.
- The locking loop pigtail aids in catheter fixation.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gage
G07340	ULT16.0-38-25-P-6S-LMCL*	16	.038	25	6	18
G07774	ULT18.0-38-25-P-6S-LMCL*	18	.038	25	6	18
G08226	ULT20.0-38-25-P-6S-LMCL*	20	.038	25	6	18

*Nonstock

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Thal-Quick Abscess Drainage Set

Used for abscess drainage.

- Multiple large, oval sideports increase drainage capacity.
- The large-bore, clear polyvinylchloride catheter features a radiopaque stripe to help visualise the catheter and its sideports during placement under fluoroscopy.
- The EchoTip® needle cannula enhances visibility under ultrasound.
- The heavy-duty wire guide offers increased support to facilitate placement of the catheter.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Dilators Fr	Needle gage/Length cm
G06336	TQAS-1200-J	12	.038	33	9	6, 10, 13	18/15
G08454	TQAS-1400-J	14	.038	41	9	6, 10, 15	18/15
G06339	TQAS-1600-J	16	.038	41	9	8, 12, 18	18/15
G06343	TQAS-1800-J	18	.038	41	9	8, 14, 20	18/15
G06345	TQAS-2000-J	20	.038	41	9	8, 14, 22	18/15
G06348	TQAS-2400-J	24	.038	41	9	12, 18, 26	18/15
G06354	TQAS-2800-J	28	.038	41	9	14, 20, 30	18/15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Gordon Large-Bore Drainage Catheter

Used for abscess drainage.


- The large-bore catheter facilitates drainage of thick fluids that will not adequately drain through smaller devices.
- Ultrathane® material enhances patient comfort.
- The catheter can be inserted by Seldinger access technique. Access components are not supplied with the 22.0 Fr coaxial catheter.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during catheter insertion.

Components

Standard catheter components




Coaxial catheter components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gage
Standard Catheter						
G12509	ULT16.0-38-40-P-6S-MPA-GLD-HC	16	.038	40	6	18
G12510	ULT18.0-38-40-P-6S-MPA-GLD-HC	18	.038	40	6	18
G12511	ULT20.0-38-40-P-6S-MPA-GLD-HC	20	.038	40	6	18
Coaxial Catheter						
G12512	ULT22.0-38-40-P-6S-MPA-GLD-HC	22	.038	40	6	-

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Ring-McLean Sump Drainage Set – Straight

Used for abscess drainage.

- The sump lumen allows the free flow of air and provides the ability to flush the catheter while the drain lumen is draining fluid.
- The straight tip configuration facilitates drainage in small cavities.
- Ultrathane® material enhances patient comfort.
- The catheter can be inserted by direct puncture with the trocar stylet that is included or by Seldinger access technique.

Components



Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gage
G07505	RMSU-12-30	12	30	13	15
G07224	RMSU-14-30	14	30	13	15
G07504	RMSU-16-30	16	30	13	15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Ring-McLean Sump Drainage Set – Curved

Used for abscess drainage.

- The sump lumen allows the free flow of air and provides the ability to flush the catheter while the drain lumen is draining fluid.
- The curved tip configuration facilitates drainage in small cavities.
- Ultrathane® material enhances patient comfort.
- The catheter can be inserted by direct puncture with the trocar styllet that is included or by Seldinger access technique.

Components



Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Trocar Styllet Needle gage
G07725	RMSU-12-30-MPB	12	30	5	15
G07726	RMSU-14-30-MPB	14	30	5	15
G07727	RMSU-16-30-MPB	16	30	5	15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Cook-Cope Ring-McLean Sump Drainage Set – Locking Pigtail

Used for abscess drainage.

- The sump lumen allows the free flow of air and provides the ability to flush the catheter while the drain lumen is draining fluid.
- The locking loop pigtail aids in catheter fixation.
- Ultrathane® material enhances patient comfort.
- The catheter can be inserted by direct puncture with the trocar stylet that is included or by Seldinger access technique.

Components



Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gage
G07768	RMSU-12-18-ACL	12	18	8	15
G07767	RMSU-14-18-ACL	14	18	8	15
G07771	RMSU-16-18-ACL	16	18	8	15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Notes

[illegible]

Notes

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Notes

[illegible]

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, is.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4