


# Physician choice, patient focus.

---




## Amplatz

URETERAL STENT SET

The technology that only a leader can provide, including Ultrathane® construction, hydrophilic coating for ease of insertion, and radiopacity for precise placement

The product quality that you have come to expect, such as kink-resistance

The radiopaque bands that indicate the location of the proximal and distal loops, enhancing visibility for accurate positioning




# Amplatz

## URETERAL STENT SET

The device, which consists of a ureteral stent and components for placement, is delivered percutaneously and is intended to establish drainage from the ureteropelvic junction to the bladder via stenting of the ureter for adult patients in whom it is desirable to place a drain that does not extend externally.

The stent is not intended to remain indwelling for more than 180 days.

### Components



| Order Number | Reference Part Number | Stent Fr | Accepts Wire Guide Diameter inch | Stent Length cm | Stent Sideports |
|--------------|-----------------------|----------|----------------------------------|-----------------|-----------------|
| G09989       | UTSSW-8.5-12-AMP-RH   | 8.5      | 0.035                            | 12              | 6               |
| G09865       | UTSSW-8.5-20-AMP-RH   | 8.5      | 0.035                            | 20              | 6               |
| G09708       | UTSSW-8.5-22-AMP-RH   | 8.5      | 0.035                            | 22              | 6               |
| G09709       | UTSSW-8.5-24-AMP-RH   | 8.5      | 0.035                            | 24              | 6               |
| G09710       | UTSSW-8.5-26-AMP-RH   | 8.5      | 0.035                            | 26              | 6               |
| G09711       | UTSSW-8.5-28-AMP-RH   | 8.5      | 0.035                            | 28              | 6               |
| G09713       | UTSSW-10.2-24-AMP-RH  | 10.2     | 0.038                            | 24              | 6               |

Please follow any step-by-step instructions carefully, including the risk, technical information, Device Name and intended use provided in the IFU at [cookmedical.eu](http://cookmedical.eu). Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.



### Customer Service

EU website: [cookmedical.eu](http://cookmedical.eu)  
USA website: [cookmedical.com](http://cookmedical.com)