

Your measure of visibility

From its highly visible tip technology to its gold sizing markers, you can expect unparalleled visibility for vessel and device measurement.

VISIBILITY:

Our iconic Beacon Tip technology allows for enhanced visibility, simplifying catheter positioning.

PRECISION:

Brilliant radiopaque marker bands ensure visibility for accurate measurement and device selection.

QUALITY:


Unique foil packaging is designed to protect the integrity of the catheter and withstand environmental factors.

Beacon[®] Tip
SIZING CATHETER



Beacon® Tip Sizing Catheter


Intended for use in angiographic procedures by physicians trained and experienced in angiographic techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed. Beacon Tip Sizing Catheters have marker bands that can be used for anatomical measurements.



Cava Vessel Sizing Catheter



Centimeter Sizing Catheter



White Vessel Sizing Catheter

Order Number	Reference Part Number	Fr	Accepts Wire Guide	Length cm	Sideports	Tip Configuration	Maximum Flow Rate mL/sec	Maximum Recommended Pressure		
			Diameter inch					kg/cm ²	psi	
Cava Vessel Sizing Catheter (Two gold bands spaced 28 mm apart)										
G31230	NR5.0-35-70-P-10S-PIG-CAVA	5.0	0.035	70	10	PIG	21	69	1000	
G31232	NR5.0-35-70-P-10S-VCF-CAVA	5.0	0.035	70	10	VCF	21	68	1000	
G31231	NR5.0-35-90-P-10S-PIG-CAVA	5.0	0.035	90	10	PIG	17	68	1000	
G31233	NR5.0-35-90-P-10S-VCF-CAVA	5.0	0.035	90	10	VCF	17	69	1000	
Centimeter Sizing Catheter (Gold bands mark 1 cm intervals along a 20 cm segment)										
G31215	NR5.0-35-70-P-10S-PIG-CSC-20	5.0	0.035	70	10	PIG	21	69	1000	
G31213	NR5.0-35-70-P-10S-0-CSC-20	5.0	0.035	70	10	Straight	21	69	1000	
G31220	NR5.0-35-70-P-10S-VCF-CSC-20	5.0	0.035	70	10	VCF	21	69	1000	
G31216	NR5.0-35-100-P-10S-PIG-CSC-20	5.0	0.035	100	10	PIG	17	69	1000	
G31214	NR5.0-35-100-P-10S-0-CSC-20	5.0	0.035	100	10	Straight	17	69	1000	
G31221	NR5.0-35-100-P-10S-VCF-CSC-20	5.0	0.035	100	10	VCF	17	69	1000	
White Vessel Sizing Catheter (Two gold bands spaced 2 cm apart)										
G31229	NR5.0-35-65-P-6S-VCF-WSC	5.0	0.035	65	6	VCF	20	69	1000	
G31224	NR5.0-35-70-P-10S-PIG-WSC	5.0	0.035	70	10	PIG	21	69	1000	
G31228	NR5.0-35-70-P-10S-0-WSC	5.0	0.035	70	10	Straight	21	69	1000	
G31226	NR5.0-35-100-P-10S-PIG-WSC	5.0	0.035	100	10	PIG	17	69	1000	

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-8.5X11