


Built for precise below-the-knee access.

Smooth hydrophilic performance, supportive shaft control, and distal visibility designed to help you navigate tortuous peripheral anatomy with confidence.



Approach[®] Hydro ST

HYDROPHILIC MICRO WIRE GUIDE

Engineered for peripheral trackability


The short taper design and hydrophilic distal segment are optimised for smooth advancement, helping maintain pushability and trackability in smaller peripheral vessels.

Visibility for precise wire control

The distal platinum microcoil helps enhance tip visualisation so physicians can confidently track wire position and make controlled movements during complex peripheral interventions.

Support where distal access demands it

A supportive stainless-steel shaft helps deliver the stability needed to smoothly advance support catheters through peripheral vessels, even in challenging below-the-knee anatomy.



Approach® Hydro ST

HYDROPHILIC MICRO WIRE GUIDE

The Approach Hydro ST Hydrophilic Micro Wire Guide facilitates delivery of devices into the peripheral vasculature.

Order Number	Reference Part Number	Wire Guide Diameter inch	Wire Guide Length cm
G52937	HMW-14-135-ST	0.014	135
G52938	HMW-14-190-ST	0.014	190
G52939	HMW-14-300-ST	0.014	300

Please follow any step-by-step instructions carefully, including the risk, technical information, Device Name and intended use provided in the IFU at cookmedical.eu. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

The following products are compatible with Approach Hydro ST.

ACCESS

Micropuncture®
Introducer Set

Micropuncture
Pedal Access Set

TARGET

CXI®
Support Catheter

Flexor® Ansel Guiding Sheath
with High-Flex Dilator

TREAT

Advance® 14LP
Low-Profile PTA Balloon
Dilatation Catheter

Customer Service

EU website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, emea.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 2161 272 3770, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com
 EDI: cookmedical.com/edi.do

Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 Email: customersupport@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 Email: cau.custserv@cookmedical.com



AI, ESC, IR, PI-A4