


Doctor choice, patient focus



PERCUTANEOUS ABSCESS AND MULTIPURPOSE DRAINAGE PRODUCTS

Cook Medical offers one of the widest assortments of drainage products available.

In this catalogue you will find the following:

- The right choice of multipurpose drainage catheters to suit each procedure
- The technology that only a leader can provide, including Ultrathane® construction, hydrophilic coating for ease of insertion, and radiopacity for precise placement
- The product qualities that you have come to expect, such as kink resistance and consistent hub sizes

This catalogue focuses on our percutaneous abscess and multipurpose drainage products.

To see our full range of percutaneous drainage products, please request these additional catalogues from your Cook representative:

- **Percutaneous Biliary Intervention Products**
- **Percutaneous Urinary Drainage Products**
- **Percutaneous Access, Targeting and Accessory Products**

Be sure to have the right drainage devices for each procedure.

The choice is yours.

Contents

Non-locking multipurpose drainage

Drainage Catheter Needle Set	4
Universal Curved Drainage Catheter	5
Multipurpose Drainage Catheter	6

Locking multipurpose drainage

Dawson-Mueller Drainage Catheter	7
Multipurpose Drainage Catheter Set – Micropuncture Set	8
Multipurpose Drainage Catheter	9

Large lumen drainage

Multipurpose Drainage Catheter – Cook-Cope Locking Loop	10
Thal-Quick Abscess Drainage Set	11
Gordon Large-Bore Drainage Catheter	12

Double lumen drainage

Ring-McLean Sump Drainage Set – Straight	13
Ring-McLean Sump Drainage Set – Curved	14
Ring-McLean Sump Drainage Set – Cook-Cope Locking Pigtail	15






Drainage Catheter Needle Set

Used for percutaneous external drainage in multiple drainage applications (e.g., nephrostomy, abscess drainage, and other abdominal cavity drainage procedures).

- The catheter features a non-locking pigtail for retention.
- The trocar stylet fits into the catheter and the stiffening cannula. This assembly facilitates access by direct puncture.
- The catheter is supplied with a pigtail straightener.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gauge
G56888	DUAN-5.0-19.5-15.0-6S-PIG	5	0.038	15	6	19.5
G56889	DUAN-5.0-19.5-25.0-6S-PIG	5	0.038	25	6	19.5
G56890	DUAN-6.0-19.5-25.0-6S-PIG	6	0.038	25	6	19.5
G56891	DUAN-7.0-19.5-25.0-8S-PIG	7	0.038	25	8	19.5
G56892	DUAN-8.3-18-25.0-8S-PIG	8.3	0.038	25	8	18

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Universal Curved Drainage Catheter

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess), either by direct puncture or Seldinger access technique.

- The curved tip configuration facilitates drainage in small cavities.
- The catheter is made of Ultrathane® material.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during insertion.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gauge
G10956	ULT6.0-35-25-P-3S-UCD-HC	6	0.035	25	3	19 UTW*
G10957	ULT7.0-38-25-P-3S-UCD-HC	7	0.038	25	3	18
G10958	ULT8.5-38-25-P-3S-UCD-HC	8.5	0.038	25	3	18
G10959	ULT10.2-38-25-P-3S-UCD-HC	10.2	0.038	25	3	18

*Ultrathin wall

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Multipurpose Drainage Catheter

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess), either by direct puncture or Seldinger access technique.

- The catheter features a non-locking pigtail for retention.
- The catheter is made of Ultrathane® material.
- The catheter can be inserted by direct puncture with the trocar stylet that is included or by Seldinger access technique.


Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gauge
G06483	ULT8.5-38-25-P-6S-MDC*	8.5	0.038	25	6	18
G06484	ULT10.2-38-25-P-6S-MDC	10.2	0.038	25	6	18
G06485	ULT12.0-38-25-P-6S-MDC	12	0.038	25	6	18
G06482	ULT14.0-38-25-P-6S-MDC	14	0.038	25	6	18

*Features Intro-Tip™

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Dawson-Mueller Drainage Catheter

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess), either by direct puncture or Seldinger access technique.

- The 10 mm pigtail facilitates drainage in small cavities.
- The catheter is made of Ultrathane® material.
- The Mac-Loc® locking mechanism firmly secures the pigtail.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during catheter insertion.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gauge
Catheter Standard Length						
G11018	ULT5.0-18-25-P-5S-CLDM-HC	5	0.018	25	5	23 UTW*
G51594	ULT6.3-35-15-P-5S-CLDM-HC	6.3	0.035	15	5	19 UTW*
G51595	ULT6.3-35-25-P-5S-CLDM-HC	6.3	0.035	25	5	19 UTW*
G11020	ULT7.0-35-25-P-5S-CLDM-HC	7	0.035	25	5	19 UTW*
G11386	ULT8.5-38-15-P-5S-CLDM-HC	8.5	0.038	15	5	18
G09705	ULT8.5-38-25-P-5S-CLDM-HC	8.5	0.038	25	5	18
G11387	ULT10.2-38-15-P-5S-CLDM-HC	10.2	0.038	15	5	18
G09706	ULT10.2-38-25-P-5S-CLDM-HC	10.2	0.038	25	5	18
G11720	ULT12.0-38-15-P-5S-CLDM-HC	12	0.038	15	5	18
G10435	ULT12.0-38-25-P-5S-CLDM-HC	12	0.038	25	5	18
G10436	ULT14.0-38-25-P-5S-CLDM-HC	14	0.038	25	5	18

*Ultrathin wall

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Multipurpose Drainage Catheter Set – Micropuncture Set

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess).

- Catheter placement is performed using micropuncture with the Neff Percutaneous Access Set that is conveniently included.
- The catheter is made of Ultrathane® material.
- The Mac-Loc® locking mechanism firmly secures the 25 mm pigtail.
- The radiopaque band identifies the proximal area of the locking loop, which contains six large oval sideports to facilitate drainage.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during catheter insertion.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Dilators Fr
G10049	CLM-8.5-RH-NPAS-NT*	8.5	0.038	25	6	9
G10050	CLM-10.2-RH-NPAS-NT	10.2	0.038	25	6	8, 11
G10261	CLM-12.0-RH-NPAS-NT	12	0.038	25	6	10, 14

*Features Intro-Tip™

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Multipurpose Drainage Catheter

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess), either by direct puncture or Seldinger access technique.

- The catheter is made of Ultrathane® material.
- The Mac-Loc® locking mechanism firmly secures the 25 mm pigtail.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during catheter insertion.
- The radiopaque band identifies the proximal area of the locking loop, which contains six large oval sideports to facilitate drainage.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gauge
G11105	ULT8.5-38-15-P-6S-CLM-RH*	8.5	0.038	15	6	18
G09501	ULT8.5-38-25-P-6S-CLM-RH*	8.5	0.038	25	6	18
G09765	ULT8.5-38-45-P-6S-CLM-RH*	8.5	0.038	45	6	18
G11210	ULT8.5-38-60-P-6S-CLM-RH* **	8.5	0.038	60	6	18
G09502	ULT10.2-38-25-P-6S-CLM-RH	10.2	0.038	25	6	18
G09766	ULT10.2-38-45-P-6S-CLM-RH	10.2	0.038	45	6	18
G28916	ULT10.2-38-60-P-6S-CLM-RH**	10.2	0.038	60	6	18
G09503	ULT12.0-38-25-P-6S-CLM-RH	12	0.038	25	6	18
G09767	ULT12.0-38-45-P-6S-CLM-RH	12	0.038	45	6	18
G09504	ULT14.0-38-25-P-6S-CLM-RH	14	0.038	25	6	18
G09768	ULT14.0-38-45-P-6S-CLM-RH**	14	0.038	45	6	18

*Features Intro-Tip™

**Nonstock

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Multipurpose Drainage Catheter – Cook-Cope Locking Loop

Used for percutaneous drainage in a variety of drainage applications (e.g., nephrostomy, biliary, and abscess), either by direct puncture or Seldinger access technique.

- The large-bore catheter features a large lumen for higher flow.
- The catheter is made of Ultrathane® material.
- The locking loop pigtail aids in catheter fixation.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gauge
G07774	ULT18.0-38-25-P-6S-LMCL*	18	0.038	25	6	18

*Nonstock

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Thal-Quick Abscess Drainage Set

Used for abscess drainage.

- Multiple large, oval sideports increase drainage capacity.
- The large-bore, clear polyvinylchloride catheter features a radiopaque stripe to help visualise the catheter and its sideports during placement under fluoroscopy.
- The EchoTip® needle cannula enhances visibility under ultrasound.
- The heavy-duty wire guide offers increased support to facilitate placement of the catheter.

Components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Dilators Fr	Needle gauge/Length cm
G06336	TQAS-1200-J	12	0.038	33	9	6, 10, 13	18/15
G08454	TQAS-1400-J	14	0.038	41	9	6, 10, 15	18/15
G06339	TQAS-1600-J	16	0.038	41	9	8, 12, 18	18/15
G06343	TQAS-1800-J	18	0.038	41	9	8, 14, 20	18/15
G06345	TQAS-2000-J	20	0.038	41	9	8, 14, 22	18/15
G06348	TQAS-2400-J	24	0.038	41	9	12, 18, 26	18/15
G06354	TQAS-2800-J	28	0.038	41	9	14, 20, 30	18/15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Gordon Large-Bore Drainage Catheter

Used for abscess drainage.


- The catheter is made of Ultrathane® material.
- The catheter can be inserted by Seldinger access technique. Access components are not supplied with the 22.0 Fr coaxial catheter.
- The hydrophilic coating becomes very lubricious when activated, which reduces friction during catheter insertion.

Components

Standard catheter components



Coaxial catheter components



Order Number	Reference Part Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gauge
Standard Catheter						
G12509	ULT16.0-38-40-P-6S-MPA-GLD-HC	16	0.038	40	6	18
G12510	ULT18.0-38-40-P-6S-MPA-GLD-HC	18	0.038	40	6	18
G12511	ULT20.0-38-40-P-6S-MPA-GLD-HC	20	0.038	40	6	18
Coaxial Catheter						
G12512	ULT22.0-38-40-P-6S-MPA-GLD-HC	22	0.038	40	6	-

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Ring-McLean Sump Drainage Set – Straight

Used for abscess drainage.

- The sump lumen allows the free flow of air and provides the ability to flush the catheter while the drain lumen is draining fluid.
- The straight tip configuration facilitates drainage in small cavities.
- The catheter is made of Ultrathane® material.
- The catheter can be inserted by direct puncture with the trocar stylet that is included or by Seldinger access technique.

Components



Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gauge
G07505	RMSU-12-30	12	30	13	15
G07224	RMSU-14-30	14	30	13	15
G07504	RMSU-16-30	16	30	13	15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Ring-McLean Sump Drainage Set – Curved

Used for abscess drainage.

- The sump lumen allows the free flow of air and provides the ability to flush the catheter while the drain lumen is draining fluid.
- The curved tip configuration facilitates drainage in small cavities.
- The catheter is made of Ultrathane® material.
- The catheter can be inserted by direct puncture with the trocar stylet that is included or by Seldinger access technique.

Components



Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gauge
G07725	RMSU-12-30-MPB	12	30	5	15
G07726	RMSU-14-30-MPB	14	30	5	15
G07727	RMSU-16-30-MPB	16	30	5	15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.




Ring-McLean Sump Drainage Set – Cook-Cope Locking Pigtail

Used for abscess drainage.

- The sump lumen allows the free flow of air and provides the ability to flush the catheter while the drain lumen is draining fluid.
- The locking loop pigtail aids in catheter fixation.
- The catheter is made of Ultrathane® material.
- The catheter can be inserted by direct puncture with the trocar stylet that is included or by Seldinger access technique.

Components



Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gauge
G07768	RMSU-12-18-ACL	12	18	8	15
G07767	RMSU-14-18-ACL	14	18	8	15
G07771	RMSU-16-18-ACL	16	18	8	15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at cookmedical.eu.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, is.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, PI-A4