

Remain in
control.



The Pin Vise Polymer is compatible for use with the Byrd Dilator Sheath Set.

Pin Vise Polymer

The Pin Vise Polymer is a chuck-like device that can be advanced over locking stylets. Tightening down the rotational mechanism at its proximal end causes the collet to close inward and grasp the outer diameter of any sheath for easier manipulation.





Pin Vise Polymer

Intended to aid the gripping of a dilator sheath by attaching to the proximal end of the sheath.

Order Number	Reference Part Number
G58341	LR-PVP001

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Byrd Dilator Sheaths Telescoping Polypropylene

Order Number	Reference Part Number	Inner Sheath ID/OD Fr	Outer Sheath ID/OD Fr	Length Inner/Outer cm
Regular				
G26562	LR-PPLBES-7.0	7.0/9.3	10.0/12.3	38/33
G19893	LR-PPLBES-8.5	8.5/10.7	11.6/13.9	38/33
G19859	LR-PPLBES-10.0	10.0/12.1	13.1/15.2	38/33
G19860	LR-PPLBES-11.5	11.5/13.6	14.1/16.3	38/33
G26514	LR-PPLBES-13.0	13.0/15.3	16.0/18.1	38/33
Extra long				
G26563	LR-PPLBES-7.0-XL	7.0/9.3	10.0/12.3	46/41
G20534	LR-PPLBES-8.5-XL	8.5/10.7	11.6/13.9	46/41
G20535	LR-PPLBES-10.0-XL	10.0/12.1	13.1/15.2	46/41
G20525	LR-PPLBES-11.5-XL	11.5/13.6	14.1/16.3	46/41
G20526	LR-PPLBES-13.0-XL	13.0/15.3	16.0/18.1	46/41

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Please see product risk information in the IFU at cookmedical.eu.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4