

Behold kontrollen.



Stiftskruepolymeren er kompatibel for bruk med Byrd dilatatorhylsesett.

Stiftskruepolymer

Stiftskruepolymeren er en chucklignende anordning som kan føres over låsestiletter. Når rotasjonsmekanismen strammes til i den proksimale enden, lukkes spennhylsen innover og griper tak i ytterdiameteren på en hylse, uansett type, for enklere håndtering.





Stiftskruepolymer

Tiltenkt å gjøre det lettere å gripe en dilatatorhylse ved å feste den til den proksimale enden av hylsen.

Bestillingsnummer Referansenummer

G58341 LR-PVP001

Enkelte produkter eller varenumre er kanskje ikke tilgjengelige i alle markeder. Kontakt din lokale Cook-representant eller kundestøtte og distribusjon for nærmere opplysninger.

Byrd teleskopisk polypropylen til dilatatorhylser

Bestillingsnummer	Referansenummer	Innvendig hylse ID/YD Fr	Utvendig hylse ID/YD Fr	Lengde innvendig/ utvendig cm
Normal				
G26562	LR-PPLBES-7.0	7,0/9,3	10,0/12,3	38/33
G19893	LR-PPLBES-8.5	8,5/10,7	11,6/13,9	38/33
G19859	LR-PPLBES-10.0	10,0/12,1	13,1/15,2	38/33
G19860	LR-PPLBES-11.5	11,5/13,6	14,1/16,3	38/33
G26514	LR-PPLBES-13.0	13,0/15,3	16,0/18,1	38/33
Ekstra lang				
G26563	LR-PPLBES-7.0-XL	7,0/9,3	10,0/12,3	46/41
G20534	LR-PPLBES-8.5-XL	8,5/10,7	11,6/13,9	46/41
G20535	LR-PPLBES-10.0-XL	10,0/12,1	13,1/15,2	46/41
G20525	LR-PPLBES-11.5-XL	11,5/13,6	14,1/16,3	46/41
G20526	LR-PPLBES-13.0-XL	13,0/15,3	16,0/18,1	46/41

Enkelte produkter eller varenumre er kanskje ikke tilgjengelige i alle markeder. Kontakt din lokale Cook-representant eller kundestøtte og distribusjon for nærmere opplysninger.

Se produktets risikoinformasjon i bruksanvisningen på cookmedical.eu.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, is.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4