





Elevate the palate. Alleviate snoring.^{1,2}



Elevo[®]

SNORING INTERVENTION SET

Effective procedure

After the lifting and stiffening action created in the soft palate, patients experienced a 30.3% reduction in snoring Visual Analog Scale at 180 days post-procedure.¹

Resorbable material

The Elevo[®] suture implant is fully resorbable, minimizing the risk of extrusions that may occur with non-resorbable implants.¹

Office-based approach


Elevoplasty[®] provides a minimally invasive, office-based procedure for patients suffering from primary snoring.



Manufactured by  Zelegant[®]

Elevo®

SNORING INTERVENTION SET



Each Elevo Snoring Intervention Set contains:

- Three preloaded Elevo Snoring Intervention Devices
- One PHIT® (Pilot Hole Initiation Tool)

Order Number	Reference Part Number	Suture Implant Nominal Length inch	Suture Implant Nominal Thickness inch	Delivery Needle Nominal Diameter inch	Quantity
G57611	EV-4000	0.76	0.023	0.079	1 set

Some products or part numbers may not be available in all markets.
Contact your local Cook Medical sales representative or Customer Service for details.

Elevo® Snoring Intervention Set

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

INDICATIONS FOR USE: The Elevo® Kit Snoring Intervention Device is intended for use in stiffening the soft palate tissue, which may reduce the severity of snoring in some individuals.

WARNINGS: Do not resterilize. Dispose of used Elevo® suturing needles in a sharps container similar to any other medically-used sharp device. Retain opened, unused Elevo® Kits only until their printed expiration date. • Store unopened Elevo® Kits in a cool place (e.g., an office that has consistent air conditioning, including on weekends) and out of direct sunlight. Storing Elevo® Kits in an office in which air conditioning is turned off on weekends may shorten shelf life and is not recommended.

PRECAUTIONS: The Elevo® Kit Snoring Intervention Device is not to be used in the treatment of obstructive sleep apnea.

ADVERSE EVENTS / COMMON TREATMENT EVENTS: Minor Bleeding • Mild Pain or Discomfort • Palate Swelling • Palpable Suture • Dysphagia

See instructions for use for full product information.

AB_50006_REV2

References

1. Elevo kit snoring intervention device. Premarket notification 510(k) K181107.
2. Indications for Use: The Elevo Kit Snoring Intervention Device is intended for use in stiffening the soft palate tissue, which may reduce the severity of snoring in some individuals.

S.I.L.E.N.C.E. clinical trial manuscript submitted for peer-review.

To learn more, visit cookmedical.com/elevo.

Elevo and Elevoplasty are registered trademarks of Zelegent, Inc.

Customer Service

EU Website: cookmedical.eu

EDI: cookmedical.eu/edi

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001702, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 972519996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Iceland: +354 800 7615, is.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4