Treat obstructive salivary gland disorders.




Illustration by Lisa Clark

SialoCath SALIVARY DUCT CATHETER

Scope-free irrigation

The catheter allows for irrigation without the use of a scope as a first-line, in-office treatment approach.

Atraumatic design

The soft tip, flexible catheter provides a gentle treatment option for both in-office and facility-based surgical procedures.

Salivary use

The catheter is specifically indicated for salivary duct irrigation.





SialoCath[™]

SALIVARY DUCT CATHETER

Used to aid in the removal and treatment of salivary duct stones and other duct conditions through flushing and/or fluid irrigation of the salivary ducts.

This flexible, single-lumen catheter can be inserted into the salivary duct either over a wire guide or through an appropriately sized introducer sheath.

Order Number	Reference Part Number	Catheter OD mm/inch	Catheter ID mm/inch	Compatible Wire Guide Diameter inch	Compatible Catheter Length cm	Tip Configuration
G35166	SDC-2.2-20	.74/.029	.53/.021	.018	15	Straight


Related products



Kolenda Salivary Access Introducer Set

Used for the introduction of diagnostic or interventional devices in the submandibular salivary duct.

Order Number	Reference Part Number	Introducer ID mm/inch	Wire Guide Diameter inch	Sheath Length cm
G24129	SDAS-6.0-18	2.2/.087	.018	5.5
G24132	SDAS-5.0-18	1.9/.074	.018	5.5



Wire guides

Used to gain and maintain access to the salivary duct.

Order Number	Reference Part Number	Wire Guide Diameter inch	Wire Guide Length inch	Tip Configuration
G02274	C-SF-15-50	.015	50	Straight
G02255	C-SF-18-40	.018	40	Straight

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details

Please see product risk information in the IFU at cookmedical.eu.

Customer Service

EU Website: cookmedical.eu EDI: cookmedical.eu/edi

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com Belgium: +32 27001702, be.orders@cookmedical.com Denmark: +45 38487607, da.orders@cookmedical.com Finland: +358 972519996, fi.orders@cookmedical.com France: +33 171230269, fr.orders@cookmedical.com Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com Iceland: +354 800 7615, IS.orders@cookmedical.com Ireland: +353 61239252, ie.orders@cookmedical.com Italy: +39 0269682853, it.orders@cookmedical.com Netherlands: +31 202013367, nl.orders@cookmedical.com Norway: +47 23162968, no.orders@cookmedical.com Spain: +34 912702691, es.orders@cookmedical.com Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335 E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com

