

# Dilate salivary duct strictures.




Illustration by Lisa Clark

## Advance<sup>®</sup>

### SALIVARY DUCT BALLOON CATHETER SET

#### High-strength

The high-strength, minimally compliant material allows for the complete dilation of the stricture.

#### Minimally invasive

The flexible, low-profile design supports minimally invasive intervention through radial dilation.

#### Multi-modal visualization

The catheter is compatible with direct endoscopic visualization or ultrasound imaging based on physician preference.





# Advance®

## SALIVARY DUCT BALLOON CATHETER SET

Used for the dilation of strictures in the submandibular and parotid salivary ducts.

The balloon is constructed of an extra-thin wall and high-strength, minimally compliant material that is intended to be used with sialendoscopic visualization in the submandibular salivary ducts or with ultrasound imaging in parotid or submandibular salivary ducts.

### Components

- Advance Salivary Duct Balloon Catheter
- Wire guide

| Order Number | Reference Part Number | Wire Guide Diameter inch | Inflated Balloon Diameter mm | Balloon Length cm | Preinflated Balloon Diameter mm |
|--------------|-----------------------|--------------------------|------------------------------|-------------------|---------------------------------|
| G34991       | SDBS-14-30-1.5-2      | .014                     | 1.5                          | 2                 | 1                               |
| G34992       | SDBS-14-30-2.0-2      | .014                     | 2.0                          | 2                 | 1                               |
| G34993       | SDBS-14-30-2.5-2      | .014                     | 2.5                          | 2                 | 1                               |
| G34994       | SDBS-14-30-3.0-2      | .014                     | 3.0                          | 2                 | 1                               |

## Related products



### Cook Sphere® Inflation Device

| Order Number | Reference Part Number | Fluid Capacity mL | Maximum Pressure atm |
|--------------|-----------------------|-------------------|----------------------|
| G34903       | CID-60-15             | 60                | 15                   |

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or the Customer Service for details.

Please see product risk information in IFU at [cookmedical.eu](http://cookmedical.eu)

### Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Iceland: +354 800 7615, [IS.orders@cookmedical.com](mailto:IS.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)  
 EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)  
**Americas:**  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)  
**Australia:**  
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



AI, ESC, IR, OHNS, PI, RH, SUR-A4