

Control and monitor balloon inflation with ease.



Cook Sphere®

INFLATION DEVICE

Easy-to-read gauge

The luminescent gauge and easy-to-read dial display pressure in both ATM and PSI.

Controlled inflation and deflation

The latch mechanism on the threaded plunger allows users to alternate between fine adjustment and rapid inflation/deflation modalities.

Enhanced ergonomics

A T-shaped handle and glove guard provide a more ergonomic grip. The spring-loaded switch allows for finer control during inflation/deflation.





Cook Sphere®

INFLATION DEVICE

Additional features

- A lower profile promotes better user control.
- Enlarged numerals on the luminescent dial and high-contrast white coloration enhance device readability and visibility in procedural lighting conditions.

Order Number	Reference Part Number	Volume mL	Maximum Pressure atm
G25571	CID-20-30-NG	20	30

Compatible Cook balloon dilatation catheters

	Wire Guide Compatibility	Balloon Inflated Diameters	Balloon Lengths	Delivery System	Catheter Length	Special Technologies
Advance 35LP	.035 inch	3–12 mm	2–20 cm	OTW	80, 135 cm	Pliaform™ texture and Starcast™ shoulders
Advance 18LP	.018 inch	2–10 mm	2–20 cm	OTW/PX	80, 135, 150, 170 cm	Rapid inflation/deflation
Advance 14LP	.014 inch	2–4 mm	2–20 cm	PX	110, 170 cm	Coiled wire shaft
Advance Micro® 14	.014 inch	1.5–3 mm	2–12 cm	OTW	50, 90, 150 cm	3 Fr compatible

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com

